



***NOW FOR
REMOVABLE & FIXED
FULL-ARCH CASES***



LOCATOR® IMPLANT SYSTEM


PRODUCT CATALOG



THE LOCATOR® IMPLANT SYSTEM

ALL-INCLUSIVE PACKAGE CONTAINS

- 1 Implant
- 1 LOCATOR Abutment
- 1 Processing Pack



2.4mm

Part	Length	Cuff
07450	10mm	2.5mm
07451	12mm	2.5mm
07452	14mm	2.5mm
07440	10mm	3mm
07441	12mm	3mm
07442	14mm	3mm
07455	10mm	4mm
07456	12mm	4mm
07457	14mm	4mm
07432	10mm	5mm
07433	12mm	5mm
07434	14mm	5mm
07381	10mm	Implant Only
07382	12mm	Implant Only
07383	14mm	Implant Only

2.9mm

Part	Length	Cuff
07460	10mm	2.5mm
07461	12mm	2.5mm
07462	14mm	2.5mm
07443	10mm	3mm
07444	12mm	3mm
07445	14mm	3mm
07465	10mm	4mm
07466	12mm	4mm
07467	14mm	4mm
07435	10mm	5mm
07436	12mm	5mm
07437	14mm	5mm
07386	10mm	Implant Only
07387	12mm	Implant Only
07388	14mm	Implant Only

3.5mm


Part	Length	Cuff
07501-02	8mm	2.5mm
07502-02	10mm	2.5mm
07503-02	12mm	2.5mm
07504-02	14mm	2.5mm
07501-03	8mm	3mm
07502-03	10mm	3mm
07503-03	12mm	3mm
07504-03	14mm	3mm
07501-04	8mm	4mm
07502-04	10mm	4mm
07503-04	12mm	4mm
07504-04	14mm	4mm
07501-05	8mm	5mm
07502-05	10mm	5mm
07503-05	12mm	5mm
07504-05	14mm	5mm
07501	8mm	Implant Only
07502	10mm	Implant Only
07503	12mm	Implant Only
07504	14mm	Implant Only

THE LOCATOR® IMPLANT SYSTEM

PROCESSING PACK



Denture
Housing



Blue Standard
Range Insert
Low Retention





Pink Standard
Range Insert
Medium Retention



Red Extended
Range Insert
Low Retention



Block-Out
Spacer



3.9mm

Part	Length	Cuff
07505-02	8mm	2.5mm
07506-02	10mm	2.5mm
07507-02	12mm	2.5mm
07508-02	14mm	2.5mm
07505-03	8mm	3mm
07506-03	10mm	3mm
07507-03	12mm	3mm
07508-03	14mm	3mm
07505-04	8mm	4mm
07506-04	10mm	4mm
07507-04	12mm	4mm
07508-04	14mm	4mm
07505-05	8mm	5mm
07506-05	10mm	5mm
07507-05	12mm	5mm
07508-05	14mm	5mm
07505	8mm	Implant Only
07506	10mm	Implant Only
07507	12mm	Implant Only
07508	14mm	Implant Only

4.4mm
































Part	Length	Cuff
07509-02	8mm	2.5mm
07510-02	10mm	2.5mm
07511-02	12mm	2.5mm
07512-02	14mm	2.5mm
07509-03	8mm	3mm
07510-03	10mm	3mm
07511-03	12mm	3mm
07512-03	14mm	3mm
07509-04	8mm	4mm
07510-04	10mm	4mm
07511-04	12mm	4mm
07512-04	14mm	4mm
07509-05	8mm	5mm
07510-05	10mm	5mm
07511-05	12mm	5mm
07512-05	14mm	5mm
07509	8mm	Implant Only
07510	10mm	Implant Only
07511	12mm	Implant Only
07512	14mm	Implant Only

4.9mm








Part	Length	Cuff
07513-02	8mm	2.5mm
07514-02	10mm	2.5mm
07515-02	12mm	2.5mm
07516-02	14mm	2.5mm
07513-03	8mm	3mm
07514-03	10mm	3mm
07515-03	12mm	3mm
07516-03	14mm	3mm
07513-04	8mm	4mm
07514-04	10mm	4mm
07515-04	12mm	4mm
07516-04	14mm	4mm
07513-05	8mm	5mm
07514-05	10mm	5mm
07515-05	12mm	5mm
07516-05	14mm	5mm
07513	8mm	Implant Only
07514	10mm	Implant Only
07515	12mm	Implant Only
07516	14mm	Implant Only


SURGICAL KIT AND INSTRUMENTATION



Part	Product	Description
07357		LOCATOR® Implant Latch Driver
07360		LOCATOR® Implant Driver, Short
07361		LOCATOR® Implant Driver, Long
07362		Torque Indicating Wrench
07463		Torque Wrench Insert
08926		LOCATOR® Square Driver Torque Wrench Insert
08340		Enhanced LOCATOR® & LOCATOR FIXED® Systems Core Tool
08536		Manual LOCATOR® & LOCATOR FIXED® Driver
07305		Drill Extender
07365		Direction Indicator
07373		LOCATOR® Implant Rotary Tissue Punch, 3.7mmD
07534		LOCATOR® Implant Rotary Tissue Punch, 5.8mmD
07537		Starter Drill, 1.5mm
07538-SB		 DLC Drill, 1.6mm
07539-SB		 DLC Drill, 2.1mm
07540-SB		 DLC Drill, 2.4mm
07541-SB		 DLC Drill, 2.6mm
07542-SB		 DLC Drill, 2.9mm
07543-SB		 DLC Drill, 3.2mm
07544-SB		 DLC Drill, 3.6mm
07545-SB		 DLC Drill, 4.1mm
07546-SB		 DLC Drill, 4.6mm

SURGICAL KIT AND INSTRUMENTATION

Part	Product	Description
07533		Narrow Ridge Drill Stops 4, 6, 8, 10, 12, 14mm
07536		Standard Ridge Drill Stops 4, 6, 8, 10, 12, 14mm
07556		Cortical Drill 3.3mm
07557		Cortical Drill 3.7mm
07558		Cortical Drill 4.3mm
07559		Cortical Drill 4.8mm
07547	Standard Surgical Tray	
07565	LOCATOR® Implant Premium Narrow Ridge Surgical Kit	
07575	LOCATOR® Implant Surgical Kit	
07576	LOCATOR® Implant Premium Surgical Kit	
09120		LOCATOR® Implant Patient Education Model

























*Items included in
Narrow Ridge Kit

GUIDED SURGICAL KIT AND INSTRUMENTATION




Part	Product	Description
07992		LOCATOR Implant Guided Implant Latch Driver
07990		LOCATOR Implant Guided Manual Implant Driver
07362		Torque Indicating Wrench
08926		LOCATOR® Square Driver Torque Wrench Insert
07463		Torque Wrench Insert
07958		LOCATOR Implant Guided Starter Drill, 1.5mmD
08340		Enhanced LOCATOR® & LOCATOR FIXED® Systems Core Tool
08536		Manual LOCATOR® & LOCATOR FIXED® Driver
08394		LOCATOR® Abutment Holder Sleeve (4 Pack)
07955		LOCATOR Implant Guided Tissue Punch, 4.9mmD
07956		LOCATOR Implant Guided Tissue Punch Starter Drill, 4.9mmD/1.5mmD
07957		LOCATOR Implant Guided Bone Profiler
07953		LOCATOR Implant Fixation Pin, Short, QTY 2
07954		LOCATOR Implant Fixation Pin, Long, QTY 2
07952		LOCATOR Implant Fixation Pin Drill
07986		LOCATOR Implant Guided Cortical Drill 3.3mm
07987		LOCATOR Implant Guided Cortical Drill 3.7mm
07988		LOCATOR Implant Guided Cortical Drill 4.3mm
07989		LOCATOR Implant Guided Cortical Drill 4.8mm
07959		● LOCATOR Implant Guided Drill, 1.6mmD x 16mmL
07960		● LOCATOR Implant Guided Drill, 1.6mmD x 20mmL
07961		● LOCATOR Implant Guided Drill, 1.6mmD x 24mmL
07962		● LOCATOR Implant Guided Drill, 2.1mmD x 16mmL
07963		● LOCATOR Implant Guided Drill, 2.1mmD x 20mmL
07964		● LOCATOR Implant Guided Drill, 2.1mmD x 24mmL


GUIDED SURGICAL KIT AND INSTRUMENTATION

Part	Product	Description
07965		● LOCATOR Implant Guided Drill, 2.4mmD x 16mmL
07966		● LOCATOR Implant Guided Drill, 2.4mmD x 20mmL
07967		● LOCATOR Implant Guided Drill, 2.4mmD x 24mmL
07968		● LOCATOR Implant Guided Drill, 2.6mmD x 16mmL
07969		● LOCATOR Implant Guided Drill, 2.6mmD x 20mmL
07970		● LOCATOR Implant Guided Drill, 2.6mmD x 24mmL
07971		● LOCATOR Implant Guided Drill, 2.9mmD x 16mmL
07972		● LOCATOR Implant Guided Drill, 2.9mmD x 20mmL
07973		● LOCATOR Implant Guided Drill, 2.9mmD x 24mmL
07974		○ LOCATOR Implant Guided Drill, 3.2mmD x 16mmL
07975		○ LOCATOR Implant Guided Drill, 3.2mmD x 20mmL
07976		○ LOCATOR Implant Guided Drill, 3.2mmD x 24mmL
07977		● LOCATOR Implant Guided Drill, 3.6mmD x 16mmL
07978		● LOCATOR Implant Guided Drill, 3.6mmD x 20mmL
07979		● LOCATOR Implant Guided Drill, 3.6mmD x 24mmL
07980		● LOCATOR Implant Guided Drill, 4.1mmD x 16mmL
07981		● LOCATOR Implant Guided Drill, 4.1mmD x 20mmL
07982		● LOCATOR Implant Guided Drill, 4.1mmD x 24mmL
07983		● LOCATOR Implant Guided Drill, 4.6mmD x 16mmL
07984		● LOCATOR Implant Guided Drill, 4.6mmD x 20mmL
07985		● LOCATOR Implant Guided Drill, 4.6mmD x 24mmL
07949	Standard Guided Surgical Tray	
07996	LOCATOR® Implant Guided Surgical Kit	
07997	LOCATOR® Implant Premium Guided Surgical Kit	
09120		LOCATOR® Implant Patient Education Model

GUIDED SURGICAL KIT AND INSTRUMENTATION




AVAILABLE IN THE
PREMIUM GUIDED SURGICAL KIT



LOCATOR® ACCESSORY COMPONENTS

MEASUREMENT COMPONENTS

Part	Product	Description
09530		Angle Measurement Guide
08517		Parallel Post (4 Pack)

PROCESSING COMPONENTS

Part	Product	Description
08505		Impression Coping (4 Pack)
08505-20		Impression Coping (20 Pack)
08530		Abutment Analog (4mm) (4 Pack)
08530-20		Abutment Analog (4mm) (20 Pack)
07339		Healing Cap, 3mm
07340		Healing Cap, 4mm
08515		Black Processing Replacement Male (4 Pack)
08515-20		Black Processing Replacement Male (20 Pack)
08514		Block-Out Spacer (20 Pack)
08569		Processing Spacer (4 Pack)
08569-20		Processing Spacer (20 Pack)


TORQUE WRENCH & DRIVERS

Part	Product	Description
09020		30Ncm Torque Wrench, Thumb, Knob and 15mm Locator Insert Driver
04391		20Ncm Torque Wrench, Thumb, Knob and 15mm Locator Insert Driver
08913		Torque Wrench Insert Driver - 23mm total length (Short) (Latch Type)
08914		Torque Wrench Insert Driver - 29mm total length (Long) (Latch Type)
08926		LOCATOR® Square Driver Torque Wrench Insert Driver (15mm)
08927		LOCATOR® Square Driver Torque Wrench Insert Driver (21mm)
04936		.050 Hex Insert Driver for 30Ncm Torque Wrench (Square Drive Connection) (21mm)
08536		Manual LOCATOR® & LOCATOR FIXED® Driver

LOCATOR® ACCESSORY COMPONENTS


STANDARD INSERTS

Dual retention to maximize stability and pivoting action that accommodates up to 20° of total divergence between two implants.








EXTENDED RANGE INSERTS

Pivoting action accommodates up to 40° of total divergence between two implants.



LOCATOR RETENTION INSERTS & DENTURE ATTACHMENT HOUSINGS

Part	Product	Description
08510		Denture Attachment Housing Assembly (4 Pack)
08510-10		Denture Attachment Housing Assembly (10 Pack)
08529		Blue: Low Retention Standard Range Inserts (4 Pack)
08529-20		Blue: Low Retention Standard Range Inserts (20 Pack)
08527		Pink: Medium Retention Standard Range Inserts (4 Pack)
08527-20		Pink: Medium Retention Standard Range Inserts (20 Pack)
08524		Clear: High Retention Standard Range Inserts (4 Pack)
08524-20		Clear: High Retention Standard Range Inserts (20 Pack)
08548		Red: Low Retention Extended Range Inserts (4 Pack)
08548-20		Red: Low Retention Extended Range Inserts (20 Pack)
08915		Orange: Medium Retention Extended Range Inserts (4 Pack)
08915-20		Orange: Medium Retention Extended Range Inserts (20 Pack)
08547		Green: High Retention Extended Range Inserts (4 Pack)
08547-20		Green: High Retention Extended Range Inserts (20 Pack)
08558		Gray: Zero Retention Extended Range Inserts (4 Pack)
08558-20		Gray: Zero Retention Extended Range Inserts (20 Pack)

LOCATOR ATTACHMENTS FOR LOCATOR IMPLANTS

Part	Description
07351	2.5mm cuff height
07318	3mm cuff height
07352	4mm cuff height
07319	5mm cuff height
07320	6mm cuff height

LOCATOR® ACCESSORY COMPONENTS

LOCATOR FIXED® RETENTION INSERTS & DENTURE ATTACHMENT HOUSINGS

Part	Product	Description
09501-2-SB		LOCATOR FIXED Insert, Four Unit (Green) (2 pack)
09501-10-SB		LOCATOR FIXED Insert, Four Unit (Green) (10 pack)
09506-2-SB		LOCATOR FIXED Insert, Anterior/Posterior (Tan) (2 pack)
09506-10-SB		LOCATOR FIXED Insert, Anterior/Posterior (Tan) (10 Pack)
09507-2-SB		LOCATOR FIXED Insert, Mid Arch (Blue) (2 pack)
09507-10-SB		LOCATOR FIXED Insert, Mid Arch (Blue) (10 pack)
09508-4-SB		LOCATOR Denture Attachment Housing Assembly, 4 pack
09508-10-SB		LOCATOR Denture Attachment Housing Assembly, 10 pack
09509-01-SB		LOCATOR FIXED Processing Package, Four Unit (Green)
09509-02-SB		LOCATOR FIXED Processing Package, Anterior/Posterior (Tan)
09509-03-SB		LOCATOR FIXED Processing Package, Mid Arch (Blue)
08340-M		Enhanced LOCATOR Systems Core Tool
09599-M		LOCATOR FIXED Seating and Removal Tool
10157-01-SB		LOCATOR FIXED Seating Tip
10156-01-SB		LOCATOR FIXED Removal Tip with Wire
10158-01-SB		LOCATOR FIXED Seating Tip Cushion
10159-01-SB		LOCATOR FIXED Removal Tip Wire Loop
10154-01-SB		LOCATOR FIXED tool 2.5mm Hex Wrench
10155-01-SB		LOCATOR FIXED Tool Tip Wrench

LOCATOR R-Tx ABUTMENTS FOR LOCATOR IMPLANTS



2.4mm

Part	Length
07381	10mm
07382	12mm
07383	14mm

2.9mm

Part	Length
07386	10mm
07387	12mm
07388	14mm

3.5mm

Part	Length
07501	8mm
07502	10mm
07503	12mm
07504	14mm

3.9mm

Part	Length
07505	8mm
07506	10mm
07507	12mm
07508	14mm



4.4mm

Part	Length
07509	8mm
07510	10mm
07511	12mm
07512	14mm

4.9mm

Part	Length
07513	8mm
07514	10mm
07515	12mm
07516	14mm

LOCATOR R-TX Abutment and Processing Pack

Part	Length
32802-02	2.5mm cuff height
32802-03	3mm cuff height
32802-04	4mm cuff height
32802-05	5mm cuff height
32802-03	6mm cuff height

THE ZEST DENTAL SOLUTIONS® 10 YEAR WARRANTY



THE ZEST DENTAL SOLUTIONS NARROW AND STANDARD RIDGE IMPLANT WARRANTY: LOCATOR® OVERDENTURE IMPLANT SYSTEM ("LODI").

Zest Dental Solutions ("Zest") is committed to providing quality products and dedicated to gathering feedback about its products. Zest actively collects and reviews the feedback of users of our products in compliance with regulatory reporting requirements and to better help us understand market expectations and validate our products' performance. The collection method Zest utilizes for such feedback is the Zest Product Experience Report ("PER") form. The PER form is to be completed from information provided by the attending clinician to share their Zest product experience.

Pursuant to the Zest implant warranty, Zest will replace covered LODI implants such qualifying covered implant is returned (as such terms are defined below). Upon a request for replacement under the Zest warranty, Zest will send a replacement product once Zest receives the returned product and completed PER form and confirms that the returned product is covered under the Zest warranty.

ONLY DIRECT ZEST CUSTOMERS MAY MAKE A WARRANTY CLAIM. In order to make a claim under the Zest warranty please return the items mentioned below in protective packaging and send these via a shipping method which enables the package to be tracked:

Printed copy of PER form completed from information provided by the attending clinician.

Explanted product(s) in sterile condition (NON-STERILE products will not qualify for replacement)

If a product failure or loss of integration has occurred, please additionally send relevant radiographs (these will not be returned unless specifically requested, please send copies)

Send shipment to: Zest Dental Solutions
ATTN: Customer Service (US customers) or Wholesale Distribution
(Distributors/OEM Partners)
2875 Loker Ave East Carlsbad, CA 92010

If using the electronic form, please include the Tracking Number for the returned product package.

Download Zest NDI Product Experience Report (PDF) AT
WWW.ZESTDENT.COM

Zest Dental Solutions ("Zest") IMPLANT WARRANTY (Valid as of May 1, 2015)

1. Warranty beneficiary and scope: Zest Dental Solutions ("Zest") hereby warrants to the direct Zest customer purchasing the from Zest ("Customer") that the Zest LOCATOR Overdenture Implant ("LODI"), when implanted according to the respective Zest LODI Technique Manual and other written instructions provided by Zest by a clinician (the "User"), will be free from any loss or lack of integration, fracture or other structural failure for the period of 10 years ("Warranty Period") from the time of treatment by the User (collectively, the "Zest Warranty"). This warranty only applies to the Customer. Third parties, particularly patients, are not covered by the Zest warranty and have no rights hereunder. Customers' sole remedy and Zest's sole liability under this Zest warranty is the replacement of the LODI implant by Zest as set forth herein. The Zest implant warranty only covers the replacement of the LODI implant and not any associated costs or expenses, including, but not limited to, chair time, laboratory fees and any other associated treatment.

2. Zest Warranty conditions: In the event that any request for warranty service for an implant is made by a Customer under the Zest warranty during the Warranty Period, Zest will replace such implant with the same or substantially equivalent product subject to the terms and conditions herein. The replacement product will be sent upon receipt of the returned product and the completed PER form and confirmation that the product is covered under the Zest warranty. To qualify for coverage under the Zest warranty, the claim must be made within the applicable Warranty Period by a Customer and all conditions below must be met. Once these requirements are satisfied, Zest will send the replacement part(s). The following REQUIRED conditions must be met and documented in order for coverage of a returned implant under

this Zest implant warranty:

- 1.** The LODI was used exclusively with all components, connections, attachments and other technology provided by Zest and not in combination with any other manufacturer's products or technology;
- 2.** The LODI is returned in sterile condition (or disinfected if delivered as such);
- 3.** The LODI was implanted by a User and inspected and maintained in full compliance with the respective Zest LODI Technique Manual valid at the time of treatment and all other Zest written instructions as well as recognized dental procedures, during and after the treatment;
- 4.** The patient had good oral hygiene which was monitored and documented by the User
- 5.** The LODI was not subjected to damage caused by misuse, misapplication, accident, trauma or any other damage caused by external factors or the User, the patient or a third party;
- 6.** A completed and signed PER was completed and submitted no later than 10 days after the complaint is made. Customer is responsible to ensure that the PER is completed and submitted from information provided by the applicable User. Details of the incident are imperative to determine whether vigilance reporting is required to regulatory authorities.

3. Limits and limitations: This Zest implant warranty is the only guarantee provided by Zest and shall apply in addition to any warranty rights conferred under any written sales agreement executed by Zest. The User remains free to claim rights against his supplier. EXCEPT AS SET FORTH HEREIN AND IN A WRITTEN SALES AGREEMENT EXECUTED BY ZEST, ZEST HEREBY DISCLAIMS ANY OTHER WARRANTIES, EXPRESS OR IMPLIED, WITH RESPECT TO LODI/SNDI OR ANY OTHER ZEST PRODUCTS, SERVICES OR INFORMATION, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTY OF FITNESS FOR PURPOSE OR MERCHANTABILITY, OR NON-INFRINGEMENT OF THIRD PARTY RIGHTS, OR ANY WARRANTY BY COURSE OF PERFORMANCE OR OTHERWISE.

ZEST SHALL NOT BE LIABLE TO THE CUSTOMER, THE USER, THE PATIENT, OR ANY THIRD PARTY, FOR LOST EARNINGS OR PROFITS, DIRECT OR INDIRECT DAMAGES AS WELL AS COLLATERAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES, DIRECTLY OR INDIRECTLY RELATED TO LODI OR ANY OTHER ZEST PRODUCTS, SERVICES OR INFORMATION.

- 4. Zest Warranty territory:** This Zest implant warranty applies worldwide to LODI implants sold by Zest, a Zest affiliated company or an official distributor of Zest on or after the validity date stated above.
- 5. Modification or termination:** Zest may modify or terminate this Zest implant warranty at any time in whole or in part. Changes to, or the termination of the Zest implant warranty, will not affect the warranty given for LODI installed prior to the date of the change or termination.

LOCATOR OVERDENTURE IMPLANT (LODI) PRODUCT EXCHANGE POLICY

Zest Dental Solutions understands that customers may need to adjust inventories of LOCATOR Overdenture Implant System (LODI) Products in order to achieve the correct mix of sizes to treat their patients. Zest Dental Solutions will waive normal restocking fees (1:1 or greater) for exchanges during the first six months following the original purchase. Packaging cannot be written on or in any way adulterated. Shipping will still be the responsibility of the customer requesting the exchange.

RETURN POLICY AND WARRANTY

ORDERING

Orders are accepted by internet, email, mail, fax, or phone. Regular business hours are 5:00am to 5:00pm PST (Pacific Standard Time). Domestic orders received by 1:00pm PST can be shipped by ground the same day. Orders received by 3:00pm PST can be shipped by express. Institutional orders require a purchase order number.

To place a new order, please choose one of the following ways:

1. Online: www.zestdent.com
2. E-mail: sales@zestdent.com
3. Phone: 1-800-262-2310
4. Fax: 1-800-487-1357

PRICING

Zest Dental Solutions makes every effort to maintain a competitive pricing structure. Pricing is subject to routine review and change without prior notice. All prices listed are in U.S. dollars.

TERMS AND BILLING

- A. Credit Card:** Payment at the time of order with VISA, AMERICAN EXPRESS, DISCOVER or MASTERCARD is necessary for customers not approved for open account billing or special financing offers.
- B. C.O.D.:** Delivery can also be made by C.O.D. freight collect.
- C. Open Account:** An open account can be established by completing a credit application and receiving approval by Zest Dental Solutions. The terms of orders purchased under open accounts are net 30 days.
- D. Export Orders:** Individual orders shipped outside the U.S. require payment in advance or a letter of credit.
- E. Past Due Accounts:** Past due balances will be subject to a 1.5% finance charge per month, amount equal to 18% per annum.

SHIPPING CHARGES

Merchandise is shipped prepaid by Zest Dental Solutions with the cost added to the invoice, or freight collect with C.O.D. fee included in the case of C.O.D. shipments. Delivery options include Fed Ex ground, 3-day, 2-day, or overnight rush service.

LIMITED WARRANTY

Zest Dental Solutions provides a limited warranty for its products, to the original purchaser, to be free from defects in workmanship and materials under normal use and service, for a period of one year from the date of purchase with the exception of LOCATOR and LOCATOR R-Tx abutments which carry a four year limited warranty (when purchased directly from Zest Dental Solutions) and a two year limited warranty when purchased from one of our partners. Zest Dental Solutions will, at its option, substitute the returned product that proves defective within the warranty period, with a similar product, free of charge.

Zest Dental Solutions continually strives to improve its products, and therefore, reserves the right to improve, modify, or discontinue products at any time without notice or incurring obligations. Purchaser assumes all risks and liability resulting from the use of Zest Dental Solutions Products, whether used separately or in combination with other products not of Zest Dental Solutions manufacture.

RETURN POLICY

Please Observe the Following Guidelines:

- A.** Authorization for returns must be received from Zest Dental Solutions prior to the return of any product. A Return Material Authorization (RMA) number will be provided for all returns.
- B.** Shipping charges must be prepaid by the customer to accept a return shipment.
- C.** Returned products are subject to a \$20 restocking fee for orders under \$400. For orders \$400 to \$1,000, a \$40 restocking fee will be applied. For orders exceeding \$1,000, please contact your Customer Service Representative for the restocking fee.
- D.** Returned merchandise will be accepted within 90 days of purchase if product is in saleable condition (in its original unopened package and not marred by any added writing or over-labeling).
- E.** Returns will not be accepted after 90 days of purchase.
- F.** Non-returnable items include:
 1. Merchandise retained beyond expiration date noted on the package.
 2. Packages with broken seals or missing parts
 3. Used, damaged, or obsolete products will not be accepted for return.
- G.** LOCATOR Implant returns will be accepted in a 1:1 exchange for a period of 6 months with no restocking fee.

VISIT OUR WEBSITE AT
WWW.ZESTDENT.COM
TO PLACE ORDERS ONLINE
24 HOURS A DAY 7 DAYS A WEEK



ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 1.855.866.LODI (5934)
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM



DISTRIBUTED BY: