



More Patients Than Ever Can Now Be Treated with LOCATOR®

INTRODUCING LOCATOR ANGLED ABUTMENT

The new LOCATOR Angled Abutment connects clinicians to a greater number of edentulous patients by offering an affordable, simple and flexible solution that meets their unique dental healthcare needs. This innovation provides 15 degrees of angle correction and makes it possible to treat more patients than ever in combination with the LOCATOR Removable or the LOCATOR FIXED restorative systems.

> DIVERGENCE CORRECTION

LOCATOR Angled Abutments correct implant angulation by 15 degrees. The LOCATOR system can restore up to an additional 20 degrees of divergence with the use of LOCATOR removable extended range or LOCATOR FIXED inserts.

> REMOVABLE OR FIXED

Compatible with LOCATOR Removable Extended Range inserts and LOCATOR FIXED Inserts.

> TREATMENT FLEXIBILITY

Angled abutments offer greater opportunity to avoid critical anatomy in implant placement and increase A-P spread to limit distal cantilevers when using LOCATOR FIXED.

> SIMPLE SOLUTION THAT REDUCES COMPLICATIONS


Properly aligned abutments not only lower risks of fractures and loosening but also ensures better retention and long-term stability. Use of this standard abutment also avoids the need for custom abutments.

> SELF-ALIGNING DESIGN

The self-aligning feature of the LOCATOR attachment allows a patient to easily seat their overdenture without the need for accurate alignment, and without causing damage to the attachment components. This self-aligning feature also increases the longevity of the LOCATOR abutment.

> TiN COATING

Features the same LOCATOR TiN coating you know and trust! LOCATOR Angled abutments allow a better path of insertion in highly angulated cases, even further protecting the abutment from wear.



REMOVABLE

*Same
LOCATOR
Abutment,
Two Solutions*



FIXED



LOCATOR Angled Abutment Ordering Information

CONNECTION	PART NUMBER	CUFF HEIGHT	REPLACEMENT SCREW	REMOVAL TOOL
Nobel Active NP Implant Direct Interactive/Simply Iconic 3.0	02415-SB	2.5mm	02414-SB	02412-SB
	02416-SB	3.5mm		
	02417-SB	4.5mm		
	02418-SB	5.5mm		
	02419-SB	6.5mm		
	02420-SB	7.5mm		
Nobel Active RP Implant Direct Interactive/Simply Iconic 3.4	02422-SB	2.5mm	02421-SB	02413-SB
	02423-SB	3.5mm		
	02424-SB	4.5mm		
	02425-SB	5.5mm		
	02426-SB	6.5mm		
	02427-SB	7.5mm		
Biohorizons Internal Hex 3.0	02462-SB	2.5mm	02461-SB	N/A
	02463-SB	3.5mm		
	02464-SB	4.5mm		
	02465-SB	5.5mm		
	02466-SB	6.5mm		
	02479-SB	7.5mm		
COMING SOON Implant Direct Legacy 3.0	02487-SB	2.5mm	02486-SB	N/A
	02488-SB	3.5mm		
	02489-SB	4.5mm		
	02490-SB	5.5mm		
	02491-SB	6.5mm		
	02492-SB	7.5mm		
ZimVie TSV 3.5 Implant Direct Legacy 3.5 Biohorizons Tapered/Internal Hex 3.5	02429-SB	2.5mm	02428-SB	N/A
	02430-SB	3.5mm		
	02431-SB	4.5mm		
	02432-SB	5.5mm		
	02433-SB	6.5mm		
	02434-SB	7.5mm		
ZimVie TSV 4.5 Implant Direct Legacy 4.5 Biohorizons Tapered/Internal Hex 4.5	02435-SB	2.5mm	02428-SB	N/A
	02436-SB	3.5mm		
	02437-SB	4.5mm		
	02438-SB	5.5mm		
	02439-SB	6.5mm		
	02440-SB	7.5mm		
Neodent GM	02442-SB	2.5mm	02441-SB	N/A
	02443-SB	3.5mm		
	02444-SB	4.5mm		
	02445-SB	5.5mm		
	02446-SB	6.5mm		
	02447-SB	7.5mm		
Straumann BLX	02398A-SB	2.5mm	02404-SB	02412-SB
	02399A-SB	3.5mm		
	02400A-SB	4.5mm		
	02401A-SB	5.5mm		
	02402A-SB	6.5mm		
	02403A-SB	7.5mm		

LOCATOR Angled Abutment Ordering Information (Continued)

CONNECTION	PART NUMBER	CUFF HEIGHT	REPLACEMENT SCREW	REMOVAL TOOL
Implant Logistics 300 Series	02450-SB	2.5mm	02449-SB	02448-SB
	02451-SB	3.5mm		
	02452-SB	4.5mm		
	02453-SB	5.5mm		
Implant Logistics 400 Series	02455-SB	2.5mm	02454-SB	02459-SB
	02456-SB	3.5mm		
	02457-SB	4.5mm		
	02458-SB	5.5mm		
Biohorizons Conical Narrow	02470-SB	2.5mm	02469-SB	02099-SB
	02471-SB	3.5mm		
	02472-SB	4.5mm		
	02480-SB	5.5mm		
	02481-SB	6.5mm		
	02482-SB	7.5mm		
Biohorizons Conical Regular	02473-SB	2.5mm	02469-SB	02099-SB
	02474-SB	3.5mm		
	02475-SB	4.5mm		
	02483-SB	5.5mm		
	02484-SB	6.5mm		
	02485-SB	7.5mm		

COMING
SOON

COMING
SOON


Required Tooling

TOOL	PART NUMBER
LOCATOR Angled Abutment Driver Set	02405-SB
LOCATOR Enhanced Core Tool	08340-M
LOCATOR Angled Abutment Removal Tool	Reference Chart Above



LOCATOR ANGLED ABUTMENT REMOVAL TOOL

LOCATOR ENHANCED CORE TOOL



LOCATOR ANGLED ABUTMENT DRIVER SET

*INCLUDES A SHORT AND LONG INSERT

CONTACT US



ZD ZEST DENTAL
SOLUTIONS®

LEARN MORE AT

WWW.ZESTDENT.COM

1.800.262.2310

SALES@ZESTDENT.COM