


LOCATOR R-Tx[®] REMOVABLE ATTACHMENT SYSTEM



U.S. RETAIL PRODUCT CATALOG


LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

System Overview and Abutment Selection Criteria

Each LOCATOR R-Tx Attachment System is delivered to you in a convenient, all-in-one double-ended vial that separately holds the Abutment and Processing Components.

Order for your favorite implant systems using this product ordering and reference guide.

LOCATOR R-Tx
Abutment




INCLUDES

- 1 - Abutment
- 1 - Denture Attachment Housing with Black Processing Insert
- 1 - Standard Zero (**Gray**) Retention Insert
- 1 - Standard Low (**Blue**) Retention Insert
- 1 - Standard Medium (**Pink**) Retention Insert
- 1 - Standard High (**Clear**) Retention Insert
- 1 - Block-Out Spacer

Standard Range Retention Inserts



Zero (**Gray**)
Retention Insert



Low (**Blue**)
Retention Insert



Medium (**Pink**)
Retention Insert



High (**Clear**)
Retention Insert

Limited Range Retention Inserts



Low (**Aqua**)
Retention Insert




Medium (**Magenta**)
Retention Insert



High (**Gold**)
Retention Insert

SELECTING THE PROPER CUFF HEIGHT FOR THE LOCATOR R-Tx ABUTMENT



1. Using a periodontal probe, measure from the top of the implant to the highest point of the gingiva.
2. Select the cuff height of the LOCATOR R-Tx Abutment that corresponds to that measurement. For example, if the measurement is 2mm, select the 2mm cuff. If the measurement is 2.5mm, select the 3mm cuff. An additional 2mm of Abutment height will extend above the gingiva to accommodate for the Denture Attachment Housing and Retention Insert.

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information

LOCATOR R-Tx ABUTMENTS ARE IMPLANT SPECIFIC

TO ORDER, YOU WILL NEED THE FOLLOWING INFORMATION:

1. Name of the implant system
2. Implant diameter
3. Connection type and platform diameter
4. Tissue measurement (at the highest point)

IMPLANT COMPANY	LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
ACE SURGICAL							
BONE LEVEL CONNECTION (Compatible with infinity Octagon Implants)							
2.7mm Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
3.3mm Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06
TISSUE LEVEL CONNECTION (Compatible with infinity Octagon Implants)							
4.8mm Regular Platform (RP)		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
INTERNAL HEX CONNECTION (Compatible with infinity Internal HEX Implants)							
3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
TRI-CHANNEL CONNECTION (Compatible with infinity TRI-CAM Implants)							
3.5mm Diameter Platform	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Diameter Platform	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Diameter Platform		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06
EXTERNAL HEX CONNECTION (Compatible with infinity EXTERNAL HEX Implants)							
3.5mm Diameter Platform		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Diameter Platform		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm Diameter Platform (Platform Switch Only)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
AVINENT							
INTERNAL HEX CONNECTION (Compatible with Biomimetic Ocean and Coral)							
3.5mm I.C (Ocean)	32500-00 (0.5mm)	32500-01	32500-02	32500-03	32500-04	32500-05	32500-06
4.1mm I.C (Ocean & Coral)	32500-00 (0.5mm)	32500-01	32500-02	32500-03	32500-04	32500-05	32500-06
EXTERNAL HEX CONNECTION (Compatible with Biomimetic Ocean and Coral)							
3.5mm E.C		32505-01	32505-02	32505-03	32505-04	32505-05	32505-06
BIOHORIZONS							
INTERNAL HEX CONNECTION (Compatible with Tapered Plus, Tapered Internal and Single-Stage)							
3.0mm		30600-01	30600-02	30600-03	30600-04	30600-05	30600-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY		LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT							
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	
BIOHORIZONS (continued)								
3.5mm	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06	
4.5mm	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06	
5.7mm		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06	
BIOMET 3i NOTE: BIOMET 3i Abutments feature a .048 hex drive								
INTERNAL HEX CERTAIN® CONNECTION (Compatible with 3i T3®, NanoTite™ and OSSEOTITE® Implants)								
3.4mm Diameter Certain Platform		30300-01	30300-02	30300-03	30300-04	30300-05	30300-06	
4.1mm Diameter Certain Platform		30301-01	30301-02	30301-03	30301-04	30301-05	30301-06	
EXTERNAL HEX CONNECTION (Compatible with 3i T3®, NanoTite™ and OSSEOTITE® Implants)								
3.4mm Diameter External Hex Platform		30305-01	30305-02	30305-03	30305-04	30305-05	30305-06	
4.1mm Diameter External Hex Platform		30306-01 (1.2mm)	30306-02	30306-03	30306-04	30306-05	30306-06	
BLUE SKY BIO								
BONE LEVEL CONNECTION (Compatible with Quattro Implants)								
2.8mm Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06	
3.3mm Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06	
TISSUE LEVEL CONNECTION (Compatible with One Stage Implants)								
4.8mm Regular Platform (RP)		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06	
6.5mm Wide Platform (WP)		30907-01	30907-02	30907-03	30907-04	30907-05	30907-06	
DENTSPLY								
CONICAL CONNECTION (Compatible with Ankylos® C/X Implants)								
2.8mm Ankylos Platform		30715-01	30715-02	30715-03	30715-04	30715-05	30715-06	
(Compatible with AstraTech Osseospeed EV Implants)								
3.6mm Diameter Platform		30700-01	30700-02	30700-03	30700-04	30700-05	30700-06	
4.2mm Diameter Platform		30701-01	30701-02	30701-03	30701-04	30701-05	30701-06	
4.8mm Diameter Platform		30702-01	30702-02	30702-03	30702-04	30702-05	30702-06	
(Compatible with AstraTech Osseospeed TX Implants)								
2.9mm Small Platform (Aqua)		30706-01	30706-02	30706-03	30706-04	30706-05	30706-06	
3.9mm Large Platform (Lilac)		30707-01	30707-02	30707-03	30707-04	30707-05	30707-06	
INTERNAL CONICAL CONNECTION								
(Compatible with OmniTaper EV, Astra Tech EV - OsseoSpeed®, and PrimeTaper EV)								
3.6mm		30700-01	30700-02	30700-03	30700-04	30700-05	30700-06	
4.2 mm		30701-01	30701-02	30701-03	30701-04	30701-05	30701-06	
4.8 mm		30702-01	30702-02	30702-03	30702-04	30702-05	30702-06	



LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY		LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)					
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
DENTSPLY (continued)							
INTERNAL HEX CONNECTION (Compatible with XiVE® Frialit-2 Implants)							
3.4mm Diameter Platform		30711-01	30711-02	30711-03	30711-04	30711-05	30711-06
3.8mm Diameter Platform		30712-01	30712-02	30712-03	30712-04	30712-05	30712-06
4.5mm Diameter Platform		30713-01	30713-02	30713-03	30713-04	30713-05	30713-06
5.5mm Diameter Platform		30714-01	30714-02	30714-03	30714-04	30714-05	30714-06
DIO							
CONICAL CONNECTION							
3.1 mm SM Narrow		31400-01	31400-02	31400-03	31400-04	31400-05	31400-06
3.5 mm SM Regular/Wide		31401-01	31401-02	31401-03	31401-04	31401-05	31401-06
2.3mm	2.3mm	31405-01	31405-02	31405-03	31405-04	31405-05	31405-06
3.35mm		31406-01	31406-02	31406-03	31406-04	31406-05	31406-06
FAIR IMPLANT							
CONICAL CONNECTION							
(Compatible with FairTwo)							
2.8mm		31900-01	31900-02	31900-03	31900-04	31900-05	31900-06
3.5mm		31901-01	31901-02	31901-03	31901-04	31901-05	31901-06
HIOSSEN/OSSTEM							
EXTERNAL HEX CONNECTION (Compatible with US)							
3.5mm Mini Platform (MP) (Yellow)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Regular Platform (RP) (Green)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
5.0mm Wide Platform Switch (WPS) (Blue)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.1mm Wide Platform (WP) (Blue)		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06
CONICAL/MORSE TAPE CONNECTION (Compatible with ET and TS)							
3.35mm Mini Regular Platform (RP) (Yellow)		31601-01	31601-02	31601-03	31601-04	31601-05	31601-06
INTERNAL HEX/MORSE TAPER CONNECTION (Compatible with ET and TS)							
2.8mm Mini Platform (MP) (Yellow)		31600-01	31600-02	31600-03	31600-04	31600-05	31600-06
INTERNAL OCTAGON/MORSE TAPER CONNECTION (Compatible with SS)							
4.8mm Regular Platform (RP) (Green)		31606-01	31606-02	31606-03	31606-04	31606-05	31606-06
IDS							
INTERNAL HEX CONNECTION (Compatible with MegaGen AnyRidge)							
3.5mm Small Platform (SP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Regular Platform (RP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY		LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)					
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
IDS (continued)							
INTERNAL HEX CONNECTION (Compatible with MegaGen AnyRidge)							
3.5mm Wide Platform (WP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
5.0mm Super Wide Platform		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
IMPLANT DIRECT							
INTERNAL HEX CONNECTION (Compatible with Legacy1, Legacy2, Legacy3 and Legacy4 Implants)							
3.0mm Diameter Platform (Blue) (Legacy 2,3,4 Only)		31500-01 (0.75mm)	31500-02	31500-03	31500-04	31500-05	31500-06
3.5mm Diameter Platform (Green)	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform (Purple)	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm Diameter Platform (Yellow)		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06
TRI-LOBE CONNECTION (Compatible with RePlant and ReActive Implants)							
3.5mm Diameter Platform (Purple)	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Diameter Platform (Yellow)	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Diameter Platform (Blue)		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06
CONICAL CONNECTION (Compatible with Interactive Implants)							
3.0mm Diameter Platform (Purple)		30506-01	30506-02	30506-03	30506-04	30506-05	30506-06
3.4mm Diameter Platform (Yellow)		30507-01	30507-02	30507-03	30507-04	30507-05	30507-06
INTERNAL OCTAGON CONNECTION (Compatible with Swish Plus and Swish Tapered Implants)							
4.8mm Diameter Platform		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06
KEYSTONE DENTAL							
CONICAL CONNECTION (Compatible with Genesis Active)							
2.9 mm KGD		31203-01	31203-02	31203-03	31203-04	31203-05	31203-06
TiLOBE CONNECTION (Compatible with Genesis and PrimaConnex)							
3.0mm Diameter (Green) Platform		31200-01	31200-02	31200-03	31200-04	31200-05	31200-06
(Compatible with PrimaConnex)							
3.6mm Diameter (Blue) Platform		31201-01	31201-02	31201-03	31201-04	31201-05	31201-06
(Compatible with PrimaConnex and TiLOBEMAXX)							
4.3mm Diameter (Purple) Platform		31202-01	31202-02	31202-03	31202-04	31202-05	31202-06
MEGAGEN/IDS							
INTERNAL HEX CONNECTION (Compatible with MegaGen AnyRidge)							
3.5mm Small Platform (SP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Regular Platform (RP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Wide Platform (WP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
5.0mm Super Wide Platform		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY		LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)					
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
MIS							
CONICAL CONNECTION (Compatible with C1 and V3 Implants)							
2.75mm Narrow Platform (NP)		32601-01	32601-02	32601-03	32601-04	32601-05	32601-06
3.15mm Standard Platform (SP)		32602-01	32602-02	32602-03	32602-04	32602-05	32602-06
INTERNAL HEX CONNECTION (Compatible with SEVEN and M4 Implants)							
3.3mm Narrow Platform (NP)		32605-01	32605-02	32605-03	32605-04	32605-05	32605-06
3.75mm Standard Platform (SP)	32606-00 (0.25mm)	32606-01	32606-02	32606-03	32606-04	32606-05	32606-06
4.5mm Wide Platform (WP)	32607-00 (0.5mm)	32607-01	32607-02	32607-03	32607-04	32607-05	32607-06
MOZO GRAU							
EXTERNAL HEX CONNECTION (Compatible with Osseous)							
4.1mm		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06
INTERNAL HEX CONNECTION (Compatible with Inhex)							
2.8mm		30706-01	30706-02	30706-03	30706-04	30706-05	30706-06
NEODENT							
CONICAL/CONE MORSE CONNECTION (Compatible with CM Implants)							
2.5mm Diameter Platform		30715-01	30715-02	30715-03	30715-04	30715-05	30715-06
CONICAL/GRAND MORSE CONNECTION (Compatible with GM Implants)							
3.0mm Diameter Platform		32700-01	32700-02	32700-03	32700-04	32700-05	32700-06
NOBEL BIOCARE							
CONICAL CONNECTION (Compatible with NobelActive®, NobelParallel™ CC and NobelReplace® CC Implants)							
3.0mm Narrow Platform (NP)		30506-01	30506-02	30506-03	30506-04	30506-05	30506-06
3.5mm Regular Platform (RP)		30507-01	30507-02	30507-03	30507-04	30507-05	30507-06
5.1mm Wide Platform (WP)		30508-01	30508-02	30508-03	30508-04	30508-05	30508-06
3.5mm Narrow Platform (NP)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Regular Platform (RP)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm Wide Platform (WP)		30512-01	30512-02	30512-03	30512-04	30512-05	
TRI-CHANNEL CONNECTION (Compatible with NobelReplace®, ReplaceSelect™ and NobelSpeedy® Replace Implants)							
3.5mm Narrow Platform (NP)	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.0mm/4.3mm Regular Platform (RP)	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Wide Platform (WP)		30502-01	30502-02	30502-03	30502-04	30502-05	

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY		LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)					
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
5.0mm Wide Platform (WP)		30502-01	30502-02	30502-03	30502-04	30502-05	
OCO BIOMEDICAL							
INTERNAL HEX CONNECTION (Compatible with Engage® Implants)							
3.0mm Diameter Platform		31500-01 (0.75mm)	31500-02	31500-03	31500-04	31500-05	31500-06
3.5mm Diameter Platform	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm Diameter Platform		30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
EXTERNAL HEX CONNECTION (Compatible with Engage® Implants)							
5.7mm		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06
OSTEOREADY/LITTLE IMPLANT COMPANY							
INTERNAL HEX CONNECTION (Compatible with Marc Nevins®, Pamela Ray™ and Performance®)							
3.5mm Diameter Platform	32401-00	32401-01	32401-02	32401-03	32401-04	32401-05	32401-06
SIC							
CONICAL CONNECTION (Compatible with SICvantage max Implants)							
2.2mm		31800-01 (1.5mm)	31800-02	31800-03	31800-04	31800-05	31800-06
INTERNAL HEX (Compatible with SICace/SICmax Implants)							
3.3mm		31805-01	31805-02	31805-03	31805-04	31805-05	31805-06
4.2mm		31806-01	31806-02	31806-03	31806-04	31806-05	31806-06
SOUTHERN IMPLANTS							
TRI-NEX CONNECTION (Compatible with Tri-Nex™ Implants)							
3.5mm Tri-Nex Platform	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Tri-Nex Platform	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Tri-Nex Platform		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06
EXTERNAL HEX CONNECTION							
4.05 Diameter Platform		31703-01	31703-02	31703-03	31703-04	31703-05	31703-06
5.0 Diameter Platform		31704-01	31704-02	31704-03	31704-04	31704-05	31704-06
6.0 Diameter Platform		31705-01	31705-02	31705-03	31705-04	31705-05	31705-06
STRAUMANN							
BONE LEVEL (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)							
3.3mm Narrow Crossfit (NC)		30401-01	30401-02	30401-03	30401-04	30401-05	30401-06
4.1mm/4.8mm Regular Crossfit (RC)		30402-01	30402-02	30402-03	30402-04	30402-05	30402-06

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM

Abutment Ordering Information (Continued)

IMPLANT COMPANY		LOCATOR R-Tx ATTACHMENT SYSTEM (Includes Abutment and Processing Components)					
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm
STRAUMANN (continued)							
(Compatible with BLX Implants)							
Regular Base (RB)		30410-01	30410-02	30410-03	30410-04	30410-05	30410-06
Wide Base (WB)		30410-01	30410-02	30410-03	30410-04	30410-05	30410-06
TISSUE LEVEL (Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)							
3.5mm Narrow Neck Crossfit (NNC)		30405-01	30405-02	30405-03	30405-04	30405-05	30405-06
4.8mm Regular Neck (RN)		30406-01	30406-02	30406-03	30406-04	30406-05	30406-06
6.5mm Wide Neck (WN)		30407-01 (1.5mm)	30407-02	30407-03	30407-04	30407-05	30407-06
SWEDEN & MARTINA							
TEXTERNAL HEX CONNECTION (Compatible with Outlink)							
3.3mm		31105-01	31105-02	31105-03	31105-04	31105-05	31105-06
4.1mm		31106-01	31106-02	31106-03	31106-04	31106-05	31106-06
INTERNAL HEX CONNECTION (Compatible with Premium/Prima)							
3.3mm		31110-01	31110-02	31110-03	31110-04	31110-05	31110-06
3.8mm		31111-01	31111-02	31111-03	31111-04	31111-05	31111-06
4.25mm		31112-01	31112-02	31112-03	31112-04	31112-05	31112-06
5.0mm		31113-01	31113-02	31113-03	31113-04	31113-05	31113-06
(Compatible with Syra)							
4.1mm		31114-01	31114-02	31114-03	31114-04	31114-05	31114-06
ZIMMER							
INTERNAL HEX CONNECTION Compatible with Tapered Screw-Vent®, Screw-Vent®, Trabecular Metal™ Implants and AdVent® 5.7mm Platform Implant)							
3.5mm	30200-00 (0.5mm)	30200-01	30200-02	30200-03	30200-04	30200-05	30200-06
4.5mm	30201-00 (0.5mm)	30201-01	30201-02	30201-03	30201-04	30201-05	30201-06
5.7mm		30202-01	30202-02	30202-03	30202-04	30202-05	30202-06
SPLINE CONNECTION (Compatible with Spline Reliance and Cylinder Implants)							
3.25mm		30205-01	30205-02	30205-03	30205-04	30205-05	30205-06
3.75mm		30206-01	30206-02	30206-03	30206-04	30206-05	30206-06
4.0mm		30206-01	30206-02	30206-03	30206-04	30206-05	30206-06
5.0mm		30207-01	30207-02	30207-03	30207-04	30207-05	30207-06





LOCATOR R-Tx[®] REMOVABLE ATTACHMENT SYSTEM

Accessory Components

PART #	DESCRIPTION	
PROCESSING COMPONENTS		
30001-01	Zero Retention Insert (Gray), (QTY. 4)	
30002-01	Low Retention Insert (Blue), (QTY. 4)	
30003-01	Medium Retention Insert (Pink), (QTY. 4)	
30004-01	High Retention Insert (Clear), (QTY. 4)	
30027-01	Low Retention Limited Range Insert (Aqua), (QTY. 4)	
30028-01	Medium Retention Limited Range Insert (Magenta), (QTY. 4)	
30029-01	High Retention Limited Range Insert (Gold), (QTY. 4)	
30012-01	Processing Insert (QTY. 4)	
08514	Block-Out Spacer (QTY. 20)	
30013-01	Denture Attachment Processing Assembly (QTY. 4)	
30014-01	3.35mm Abutment Analog (QTY. 4)	
30015-01	4mm Abutment Analog (QTY. 4)	
30016-01	5mm Abutment Analog (QTY. 4)	
30017-01	Impression Coping (QTY. 4)	
30018-01	Processing Spacer (QTY. 4)	
30030-01	Mixed Retention Limited Range Processing Pack - 2 Pack	
30031-01	Limited Range Retention Inserts, Processing Pack - 1 Pack	

LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM


Accessory Components (Continued)

TOOLS		
30021-01	LOCATOR R-Tx Retention Insert Tool (Dual Sided)	
08007	.050"/1.25mm Hex Driver with thumb knob (Shaft 6.4mm)	
08008	.050"/1.25mm Hex, Latchlock Driver (Shaft 12.7mm)	
04936	.050"/1.25mm Hex, Torque Wrench Insert (Shaft 6.4mm)	
PART #	DESCRIPTION	
PATIENT EDUCATION		
09665	LOCATOR R-Tx Patient Education Model (Straight)	
30053-01	LOCATOR R-Tx 4x Demonstration Model	
09677	LOCATOR R-Tx Patient Education Model (Angled)	
SCAN BODIES		
30060-01-2	LOCATOR Scan Body (2 Pack)	
30060-01-4	LOCATOR Scan Body (4 Pack)	
30060-01-10	LOCATOR Scan Body (10 Pack)	



09665


LOCATOR R-Tx[®]
REMOVABLE ATTACHMENT SYSTEM



VISIT **ZESTDENT.COM**
24 HOURS A DAY | 7 DAYS A WEEK

TO PLACE ORDERS ONLINE OR TO FIND
THE LATEST SUPPORT TOOLS FOR
LOCATOR R-Tx[®]

ZD ZEST DENTAL
SOLUTIONS[®]

ZEST DENTAL SOLUTIONS
2875 LOKER AVENUE EAST
CARLSBAD, CA 92010
USA

TEL: 800.262.2310
FAX: 760.743.7975
EMAIL: SALES@ZESTDENT.COM
WWW.ZESTDENT.COM

DISTRIBUTED BY: