



# 2024 PRODUCT CATALOG

ATTACHMENTS  
IMPLANTS  
MATERIALS  
EQUIPMENT

**LOCATOR  
FIXED Now  
Available!**  
See page 41 to  
learn more.



# TABLE OF CONTENTS


<b>LOCATOR® ROOT</b> .....	5
ABUTMENTS .....	6
PROCESSING COMPONENTS .....	7
TOOLS AND DRILLS .....	7
<b>LOCATOR® IMPLANTS</b> .....	9
SURGICAL INSTRUMENTS AND ACCESSORIES .....	14
<b>LOCATOR®</b> .....	18
ABUTMENTS .....	19
PROCESSING COMPONENTS .....	19
TOOLS .....	21
MULTI-UNIT ABUTMENTS .....	22
PATIENT EDUCATION MODELS .....	22
LOCATOR CONNECTIONS BY SYSTEM .....	23
<b>LOCATOR® BAR</b> .....	39
ATTACHMENT .....	39
MULTI-UNIT ABUTMENTS (BAR-SPLINTED APPLICATION) .....	40
MALE PROCESSING PACKAGES .....	40
ATTACHMENT TOOLS .....	40
<b>LOCATOR FIXED®</b> .....	41
STARTER AND COMPONENT KITS .....	42
INSERTS AND HOUSINGS .....	43
TOOLING & OTHER .....	44
PATIENT EDUCATION .....	44
<b>REGAL HYGIENE</b> .....	45
REGAL PROTECT DENTURE & FIXED-PROSTHETIC PROTECTIVE COATING .....	46
REGAL GO ALL-IN-ONE ORAL CARE BRUSH .....	46
<b>LOCATOR R-Tx®</b> .....	47
REMOVABLE ATTACHMENT SYSTEM .....	48
PROCESSING COMPONENTS .....	48
TOOLS .....	49
PATIENT EDUCATION MODELS .....	49
LOCATOR R-Tx CONNECTIONS BY SYSTEM .....	50

■ LOCATOR CHAIRSIDE® .....	58
IMPRESSION MATERIALS.....	58
BITE REGISTRATION MATERIAL .....	58
ATTACHMENT PROCESSING MATERIAL.....	58
TOOLS AND PREP AND POLISH KIT .....	59
SOFT RELINE MATERIAL KIT.....	59
BLOCK OUT MATERIAL.....	59
DENTURE REMOVAL TOOL.....	59
■ DENTAL MATERIALS & EQUIPMENT .....	60
CARIES DETECTION & DESENSITIZING.....	62
AIR INHIBITORS, BONDING AGENTS AND ETCHANTS .....	62
COMPOSITES.....	64
CEMENT.....	72
IMPRESSION MATERIAL .....	72
MATRIX SYSTEMS.....	75
STAIN REMOVAL .....	77
POSTS.....	77
PROVISIONAL COMPOSITES.....	78
IVERI WHITENING .....	80
AIR ABRASION .....	82
SANDBLASTING .....	83
ABRASIVES .....	85
DUST CONFINEMENT .....	86
RETURN POLICY AND WARRANTY .....	88

Acct Name: \_\_\_\_\_

Acct Number: \_\_\_\_\_

# MAKE IT A LOCATOR<sup>®</sup> LIFE



FROM PRE-EDENTULOUS TO FULL EDENTULISM, YOU CAN

## LIVE WITH CONFIDENCE

Discover the complete spectrum of full-arch implant solutions from Zest. See how a LOCATOR Life promotes ongoing relationships with your patients by giving them freedom to live a life with more pleasure, satisfaction and confidence.


Offer your patients outstanding aesthetics and clinically proven performance for every stage of their journey with LOCATOR systems solutions today.

## FROM REMOVABLE TO FIXED

# Attachment Systems **LOCATOR®** **ROOT**

**The gentle and easily navigable transition to dental implants.**

Leverages the natural tooth root for the retention of overdentures or removable partial dentures and offers patients all the benefits of LOCATOR without the time or cost of implants.



LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## LOCATOR® ROOT ABUTMENTS

PART#	DESCRIPTION
08011	LOCATOR Root Starter Kit (with Drills) (Includes: (4) 0° Roots, (4) Inserts Processing Packages (Denture Cap with Black Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue) Pilot Drill, Spot-Face Drill and (2) Parallel Posts)
08010	LOCATOR Root Refill Kit (No Drills) (Includes: (4) 0° Roots, (4) Inserts Processing Packages (Denture Cap with Black Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue) and (2) Parallel Posts)
08071	LOCATOR Root Mixed Kit (with Drills) (Includes: (2) 0° Roots, (1) 10° Roots, (1) 20° Roots, (4) Inserts Processing Packages (Denture Cap with Black Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue) Pilot Drill, Spot-Face Drill and (2) Parallel Posts)
08070	LOCATOR Root Mixed Refill Kit (No Drills) (Includes: (2) 0° Roots, (1) 10° Roots, (1) 20° Roots, (4) Inserts Processing Packages (Denture Cap with Black Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue) and (2) Parallel Posts)
08051	LOCATOR Root Coping Kit (with Drills) (Includes: (4) 0° Cast-To Roots, (4) Inserts Processing Packages (Denture Cap with Black Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue) Pilot Drill, Spot-Face Drill and (2) Parallel Posts)
08050	LOCATOR Root Coping Refill Kit (No Drills) (Includes: (4) 0° Cast-To Roots, (4) Inserts Processing Packages (Denture Cap with Black Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue) and (2) Parallel Posts)
01747-2	ERA Root Attachment - 2 Pack
01747-10	ERA Root Attachment - 10 Pack
08520	LOCATOR 0° 1 Pack
08520-2	LOCATOR Root 0° - 2 Pack
08520-10	LOCATOR Root 0° - 10 Pack
08521	LOCATOR 10° 1 Pack
08521-2	LOCATOR Root 10° - 2 Pack
08521-10	LOCATOR Root 10° - 10 Pack
08522	LOCATOR 20° 1 Pack

## LOCATOR® ROOT ABUTMENTS (CONTINUED)

PART#	DESCRIPTION
08522-2	LOCATOR Root 20° - 2 Pack
08522-10	LOCATOR Root 20° - 10 Pack
08528-2	LOCATOR Root 0° Cast-To Roots - 2 Pack
08528-10	LOCATOR Root 0° Cast-To Roots - 10 Pack
08522-10	LOCATOR Root 20° - 10 Pack
08528-2	LOCATOR Root 0° Cast-To Roots - 2 Pack
08528-10	LOCATOR Root 0° Cast-To Roots - 10 Pack

## LOCATOR® ROOT PROCESSING COMPONENTS

PART#	DESCRIPTION
08519-2	Inserts Processing Package - 2 Pack (Includes: Denture Cap with Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue)
08519-10	Inserts Processing Package - 10 Pack (Includes: Denture Cap with Processing Inserts, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue)

## LOCATOR® ROOT TOOLS AND DRILLS

PART#	DESCRIPTION
08340	Enhanced LOCATOR Core Tool
08393	LOCATOR Core Tool
08924	Pilot Drill
08922	Spot Face Diamond Bur

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT





# Implant Systems **LOCATOR®** **IMPLANTS**

**Economical, ridge-specific  
FIXED or Removable full-arch  
implant solution that gives  
partially or fully edentulous  
patients simplicity and  
higher quality of life.**






FIXED or Removable full-arch  
implant solution that gives  
partially or fully edentulous  
patients simplicity and higher  
quality of life.



+



PROCESSING  
PACK

		
Pink Standard Range Insert Medium Retention	Red Extended Range Insert Low Retention	Blue Standard Range Insert Low Retention
		
Denture Cap	Block-Out Spacer	

## LOCATOR® IMPLANT SYSTEM

### LOCATOR® IMPLANTS

PART#	DESCRIPTION
07450	LOCATOR Implant, <b>2.4mmD</b> x <b>10mmL</b> , <b>2.5mm</b> Cuff
07440	LOCATOR Implant, <b>2.4mmD</b> x <b>10mmL</b> , <b>3mm</b> Cuff
07455	LOCATOR Implant, <b>2.4mmD</b> x <b>10mmL</b> , <b>4mm</b> Cuff
07432	LOCATOR Implant, <b>2.4mmD</b> x <b>10mmL</b> , <b>5mm</b> Cuff
07381	LOCATOR Implant, <b>2.4mmD</b> x <b>10mmL</b> , Implant Only
07451	LOCATOR Implant, <b>2.4mmD</b> x <b>12mmL</b> , <b>2.5mm</b> Cuff
07441	LOCATOR Implant, <b>2.4mmD</b> x <b>12mmL</b> , <b>3mm</b> Cuff
07456	LOCATOR Implant, <b>2.4mmD</b> x <b>12mmL</b> , <b>4mm</b> Cuff
07433	LOCATOR Implant, <b>2.4mmD</b> x <b>12mmL</b> , <b>5mm</b> Cuff
07382	LOCATOR Implant, <b>2.4mmD</b> x <b>12mmL</b> , Implant Only
07452	LOCATOR Implant, <b>2.4mmD</b> x <b>14mmL</b> , <b>2.5mm</b> Cuff
07442	LOCATOR Implant, <b>2.4mmD</b> x <b>14mmL</b> , <b>3mm</b> Cuff
07457	LOCATOR Implant, <b>2.4mmD</b> x <b>14mmL</b> , <b>4mm</b> Cuff
07434	LOCATOR Implant, <b>2.4mmD</b> x <b>14mmL</b> , <b>5mm</b> Cuff
07383	LOCATOR Implant, <b>2.4mmD</b> x <b>14mmL</b> , Implant Only
07460	LOCATOR Implant, <b>2.9mmD</b> x <b>10mmL</b> , <b>2.5mm</b> Cuff
07443	LOCATOR Implant, <b>2.9mmD</b> x <b>10mmL</b> , <b>3mm</b> Cuff
07465	LOCATOR Implant, <b>2.9mmD</b> x <b>10mmL</b> , <b>4mm</b> Cuff
07435	LOCATOR Implant, <b>2.9mmD</b> x <b>10mmL</b> , <b>5mm</b> Cuff
07386	LOCATOR Implant, <b>2.9mmD</b> x <b>10mmL</b> , Implant Only
07461	LOCATOR Implant, <b>2.9mmD</b> x <b>12mmL</b> , <b>2.5mm</b> Cuff
07444	LOCATOR Implant, <b>2.9mmD</b> x <b>12mmL</b> , <b>3mm</b> Cuff
07466	LOCATOR Implant, <b>2.9mmD</b> x <b>12mmL</b> , <b>4mm</b> Cuff
07436	LOCATOR Implant, <b>2.9mmD</b> x <b>12mmL</b> , <b>5mm</b> Cuff
07387	LOCATOR Implant, <b>2.9mmD</b> x <b>12mmL</b> , Implant Only
07462	LOCATOR Implant, <b>2.9mmD</b> x <b>14mmL</b> , <b>2.5mm</b> Cuff
07445	LOCATOR Implant, <b>2.9mmD</b> x <b>14mmL</b> , <b>3mm</b> Cuff

**LOCATOR® IMPLANTS (CONTINUED)**

PART#	DESCRIPTION
07467	LOCATOR Implant, <b>2.9mmD</b> x <b>14mmL</b> , <b>4mm</b> Cuff
07437	LOCATOR Implant, <b>2.9mmD</b> x <b>14mmL</b> , <b>5mm</b> Cuff
07388	LOCATOR Implant, <b>2.9mmD</b> x <b>14mmL</b> , Implant Only
07501-02	LOCATOR Implant, <b>3.5mmD</b> x <b>8mmL</b> , <b>2.5mm</b> Cuff
07501-03	LOCATOR Implant, <b>3.5mmD</b> x <b>8mmL</b> , <b>3mm</b> Cuff
07501-04	LOCATOR Implant, <b>3.5mmD</b> x <b>8mmL</b> , <b>4mm</b> Cuff
07501-05	LOCATOR Implant, <b>3.5mmD</b> x <b>8mmL</b> , <b>5mm</b> Cuff
07501	LOCATOR Implant, <b>3.5mmD</b> x <b>8mmL</b> , Implant Only
07502-02	LOCATOR Implant, <b>3.5mmD</b> x <b>10mmL</b> , <b>2.5mm</b> Cuff
07502-03	LOCATOR Implant, <b>3.5mmD</b> x <b>10mmL</b> , <b>3mm</b> Cuff
07502-04	LOCATOR Implant, <b>3.5mmD</b> x <b>10mmL</b> , <b>4mm</b> Cuff
07502-05	LOCATOR Implant, <b>3.5mmD</b> x <b>10mmL</b> , <b>5mm</b> Cuff
07502	LOCATOR Implant, <b>3.5mmD</b> x <b>10mmL</b> , Implant Only
07503-02	LOCATOR Implant, <b>3.5mmD</b> x <b>12mmL</b> , <b>2.5mm</b> Cuff
07503-03	LOCATOR Implant, <b>3.5mmD</b> x <b>12mmL</b> , <b>3mm</b> Cuff
07503-04	LOCATOR Implant, <b>3.5mmD</b> x <b>12mmL</b> , <b>4mm</b> Cuff
07503-05	LOCATOR Implant, <b>3.5mmD</b> x <b>12mmL</b> , <b>5mm</b> Cuff
07503	LOCATOR Implant, <b>3.5mmD</b> x <b>12mmL</b> , Implant Only
07504-02	LOCATOR Implant, <b>3.5mmD</b> x <b>14mmL</b> , <b>2.5mm</b> Cuff
07504-03	LOCATOR Implant, <b>3.5mmD</b> x <b>14mmL</b> , <b>3mm</b> Cuff
07504-04	LOCATOR Implant, <b>3.5mmD</b> x <b>14mmL</b> , <b>4mm</b> Cuff
07504-05	LOCATOR Implant, <b>3.5mmD</b> x <b>14mmL</b> , <b>5mm</b> Cuff
07504	LOCATOR Implant, <b>3.5mmD</b> x <b>14mmL</b> , Implant Only
07505-02	LOCATOR Implant, <b>3.9mmD</b> x <b>8mmL</b> , <b>2.5mm</b> Cuff
07505-03	LOCATOR Implant, <b>3.9mmD</b> x <b>8mmL</b> , <b>3mm</b> Cuff
07505-04	LOCATOR Implant, <b>3.9mmD</b> x <b>8mmL</b> , <b>4mm</b> Cuff
07505-05	LOCATOR Implant, <b>3.9mmD</b> x <b>8mmL</b> , <b>5mm</b> Cuff

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-TX
LOCATOR CHAIRSIDE
MATERIALS & EQUIPMENT

## LOCATOR® IMPLANTS (CONTINUED)

PART#	DESCRIPTION
07505	LOCATOR Implant, <b>3.9mmD</b> x <b>8mmL</b> , Implant Only
07506-02	LOCATOR Implant, <b>3.9mmD</b> x <b>10mmL</b> , <b>2.5mm</b> Cuff
07506-03	LOCATOR Implant, <b>3.9mmD</b> x <b>10mmL</b> , <b>3mm</b> Cuff
07506-04	LOCATOR Implant, <b>3.9mmD</b> x <b>10mmL</b> , <b>4mm</b> Cuff
07506-05	LOCATOR Implant, <b>3.9mmD</b> x <b>10mmL</b> , <b>5mm</b> Cuff
07506	LOCATOR Implant, <b>3.9mmD</b> x <b>10mmL</b> , Implant Only
07507-02	LOCATOR Implant, <b>3.9mmD</b> x <b>12mmL</b> , <b>2.5mm</b> Cuff
07507-03	LOCATOR Implant, <b>3.9mmD</b> x <b>12mmL</b> , <b>3mm</b> Cuff
07507-04	LOCATOR Implant, <b>3.9mmD</b> x <b>12mmL</b> , <b>4mm</b> Cuff
07507-05	LOCATOR Implant, <b>3.9mmD</b> x <b>12mmL</b> , <b>5mm</b> Cuff
07507	LOCATOR Implant, <b>3.9mmD</b> x <b>12mmL</b> , Implant Only
07508-02	LOCATOR Implant, <b>3.9mmD</b> x <b>14mmL</b> , <b>2.5mm</b> Cuff
07508-03	LOCATOR Implant, <b>3.9mmD</b> x <b>14mmL</b> , <b>3mm</b> Cuff
07508-04	LOCATOR Implant, <b>3.9mmD</b> x <b>14mmL</b> , <b>4mm</b> Cuff
07508-05	LOCATOR Implant, <b>3.9mmD</b> x <b>14mmL</b> , <b>5mm</b> Cuff
07508	LOCATOR Implant, <b>3.9mmD</b> x <b>14mmL</b> , Implant Only
07509-02	LOCATOR Implant, <b>4.4mmD</b> x <b>8mmL</b> , <b>2.5mm</b> Cuff
07509-03	LOCATOR Implant, <b>4.4mmD</b> x <b>8mmL</b> , <b>3mm</b> Cuff
07509-04	LOCATOR Implant, <b>4.4mmD</b> x <b>8mmL</b> , <b>4mm</b> Cuff
07509-05	LOCATOR Implant, <b>4.4mmD</b> x <b>8mmL</b> , <b>5mm</b> Cuff
07509	LOCATOR Implant, <b>4.4mmD</b> x <b>8mmL</b> , Implant Only
07510-02	LOCATOR Implant, <b>4.4mmD</b> x <b>10mmL</b> , <b>2.5mm</b> Cuff
07510-03	LOCATOR Implant, <b>4.4mmD</b> x <b>10mmL</b> , <b>3mm</b> Cuff
07510-04	LOCATOR Implant, <b>4.4mmD</b> x <b>10mmL</b> , <b>4mm</b> Cuff
07510-05	LOCATOR Implant, <b>4.4mmD</b> x <b>10mmL</b> , <b>5mm</b> Cuff
07510	LOCATOR Implant, <b>4.4mmD</b> x <b>10mmL</b> , Implant Only
07511-02	LOCATOR Implant, <b>4.4mmD</b> x <b>12mmL</b> , <b>2.5mm</b> Cuff

## LOCATOR® IMPLANTS (CONTINUED)

PART#	DESCRIPTION
07511-03	LOCATOR Implant, <b>4.4mmD</b> x <b>12mmL</b> , <b>3mm</b> Cuff
07511-04	LOCATOR Implant, <b>4.4mmD</b> x <b>12mmL</b> , <b>4mm</b> Cuff
07511-05	LOCATOR Implant, <b>4.4mmD</b> x <b>12mmL</b> , <b>5mm</b> Cuff
07511	LOCATOR Implant, <b>4.4mmD</b> x <b>12mmL</b> , Implant Only
07512-02	LOCATOR Implant, <b>4.4mmD</b> x <b>14mmL</b> , <b>2.5mm</b> Cuff
07512-03	LOCATOR Implant, <b>4.4mmD</b> x <b>14mmL</b> , <b>3mm</b> Cuff
07512-04	LOCATOR Implant, <b>4.4mmD</b> x <b>14mmL</b> , <b>4mm</b> Cuff
07512-05	LOCATOR Implant, <b>4.4mmD</b> x <b>14mmL</b> , <b>5mm</b> Cuff
07512	LOCATOR Implant, <b>4.4mmD</b> x <b>14mmL</b> , Implant Only
07513-02	LOCATOR Implant, <b>4.9mmD</b> x <b>8mmL</b> , <b>2.5mm</b> Cuff
07513-03	LOCATOR Implant, <b>4.9mmD</b> x <b>8mmL</b> , <b>3mm</b> Cuff
07513-04	LOCATOR Implant, <b>4.9mmD</b> x <b>8mmL</b> , <b>4mm</b> Cuff
07513-05	LOCATOR Implant, <b>4.9mmD</b> x <b>8mmL</b> , <b>5mm</b> Cuff
07513	LOCATOR Implant, <b>4.9mmD</b> x <b>8mmL</b> , Implant Only
07514-02	LOCATOR Implant, <b>4.9mmD</b> x <b>10mmL</b> , <b>2.5mm</b> Cuff
07514-03	LOCATOR Implant, <b>4.9mmD</b> x <b>10mmL</b> , <b>3mm</b> Cuff
07514-04	LOCATOR Implant, <b>4.9mmD</b> x <b>10mmL</b> , <b>4mm</b> Cuff
07514-05	LOCATOR Implant, <b>4.9mmD</b> x <b>10mmL</b> , <b>5mm</b> Cuff
07514	LOCATOR Implant, <b>4.9mmD</b> x <b>10mmL</b> , Implant Only
07515-02	LOCATOR Implant, <b>4.9mmD</b> x <b>12mmL</b> , <b>2.5mm</b> Cuff
07515-03	LOCATOR Implant, <b>4.9mmD</b> x <b>12mmL</b> , <b>3mm</b> Cuff
07515-04	LOCATOR Implant, <b>4.9mmD</b> x <b>12mmL</b> , <b>4mm</b> Cuff
07515-05	LOCATOR Implant, <b>4.9mmD</b> x <b>12mmL</b> , <b>5mm</b> Cuff
07515	LOCATOR Implant, <b>4.9mmD</b> x <b>12mmL</b> , Implant Only
07516-02	LOCATOR Implant, <b>4.9mmD</b> x <b>14mmL</b> , <b>2.5mm</b> Cuff
07516-03	LOCATOR Implant, <b>4.9mmD</b> x <b>14mmL</b> , <b>3mm</b> Cuff
07516-04	LOCATOR Implant, <b>4.9mmD</b> x <b>14mmL</b> , <b>4mm</b> Cuff

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-TX
LOCATOR CHAIRSIDE
MATERIALS & EQUIPMENT

## LOCATOR® IMPLANTS (CONTINUED)

PART#	DESCRIPTION
07516-05	LOCATOR Implant, <b>4.9mmD</b> x <b>14mmL</b> , <b>5mm</b> Cuff
07516	LOCATOR Implant, <b>4.9mmD</b> x <b>14mmL</b> , Implant Only

## LOCATOR® ABUTMENTS FOR LOCATOR® IMPLANTS

PART#	DESCRIPTION
07351	LOCATOR Attachment for LOCATOR Implants, <b>2.5mm</b> Cuff
07318	LOCATOR Attachment for LOCATOR Implants, <b>3mm</b> Cuff
07352	LOCATOR Attachment for LOCATOR Implants, <b>4mm</b> Cuff
07319	LOCATOR Attachment for LOCATOR Implants, <b>5mm</b> Cuff
07320	LOCATOR Attachment for LOCATOR Implants, <b>6mm</b> Cuff

## LOCATOR R-TX® ABUTMENTS FOR LOCATOR® IMPLANTS

PART#	DESCRIPTION
32802-02	LOCATOR R-Tx system, LOCATOR® Implant, <b>2.5mm</b>
32802-03	LOCATOR R-Tx system, LOCATOR® Implant, <b>3mm</b>
32802-04	LOCATOR R-Tx system, LOCATOR® Implant, <b>4mm</b>
32802-05	LOCATOR R-Tx system, LOCATOR® Implant, <b>5mm</b>
32802-06	LOCATOR R-Tx system, LOCATOR® Implant, <b>6mm</b>

Go to page 42 for our compatible LOCATOR FIXED offering.

## SURGICAL INSTRUMENTS AND ACCESSORIES

PART#	DESCRIPTION
07357	LOCATOR Implant Latch Driver
07360	LOCATOR Implant Driver, Short

**SURGICAL INSTRUMENTS AND ACCESSORIES (CONTINUED)**

PART#	DESCRIPTION
07361	LOCATOR Implant, Long
07362	Torque Indicating Wrench
07463	Toque Wrench Insert
08926	LOCATOR Square Driver Toque Wrench Insert
08340	Enhanced LOCATOR Core Tool
08393	LOCATOR Core Tool
08536	Manual LOCATOR & LOCATOR FIXED Driver
08390	LOCATOR Core Tool Abutment Driver (Gold End)
08397	LOCATOR Core Tool, Insert Removal (Tip End)
08394	LOCATOR Abutment Holder Sleeve (4 Pack)
07374	Drill Extender
07365	Direction Indicator
07373	LOCATOR Implant Rotary Tissue Punch, 3.7mmD
07534	LOCATOR Implant Rotary Tissue Punch, 5.8mmD
07537	Starter Drill, 1.5mm
07538-SB	DLC Drill, 1.6mm (Green)
07539-SB	DLC Drill, 2.1mm (Cyan)
07540-SB	DLC Drill, 2.4mm (Purple)
07541-SB	DLC Drill, 2.6mm (Pink)
07542-SB	DLC Drill, 2.9mm (Yellow)
07543-SB	DLC Drill, 3.2mm (White)
07544-SB	DLC Drill, 3.6mm (Grey)
07545-SB	DLC Drill, 4.1mm (Black)
07546-SB	DLC Drill, 4.6mm (Navy)
07533	Narrow Ridge Drill Stops 4, 6, 8, 10, 12, 14mm
07536	Standard Ridge Drill Stops 4, 6, 8, 10, 12, 14mm

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-TX
LOCATOR CHAIRSIDE
MATERIALS & EQUIPMENT

## SURGICAL INSTRUMENTS AND ACCESSORIES (CONTINUED)

PART#	DESCRIPTION
07556	Cortical Drill 3.3mm (White)
07557	Cortical Drill 3.7mm (Grey)
07558	Cortical Drill 4.3mm (Black)
07559	Cortical Drill 4.8mm (Navy)
07547	Standard Surgical Tray
07997	LOCATOR Implant Premium Guided Surgical Kit
07996	LOCATOR Implant Standard Guided Surgical Kit
07575	LOCATOR Implant Surgical Kit
07576	LOCATOR Implant Premium Surgical Kit
09120	LOCATOR Implant Patient Education Model
09557	LOCATOR Overdenture 5x Implant Demonstration Model
07950-SB	LOCATOR Implant Guided Sleeve
07951-SB	LOCATOR Implant Fixation Pin Sleeve
07952-SB	LOCATOR Implant Fixation Pin Drill
07953-SB	LOCATOR Implant Fixation Pin, Short, QTY 2
07954-SB	LOCATOR Implant Fixation Pin, Long, QTY 2
07955-SB	LOCATOR Implant Guided Tissue Punch, 4.9mmD
07956-SB	LOCATOR Implant Guided Tissue Punch Starter Drill, 4.9mmD/1.5mmD
07957-SB	LOCATOR Implant Guided Bone Profiler
07958-SB	LOCATOR Implant Guided Starter Drill, 1.5mmD
07959-SB	LOCATOR Implant Guided Drill, 1.6mmD x 16mmL
07960-SB	LOCATOR Implant Guided Drill, 1.6mmD x 20mmL
07961-SB	LOCATOR Implant Guided Drill, 1.6mmD x 24mmL
07962-SB	LOCATOR Implant Guided Drill, 2.1mmD x 16mmL
07963-SB	LOCATOR Implant Guided Drill, 2.1mmD x 20mmL
07964-SB	LOCATOR Implant Guided Drill, 2.1mmD x 24mmL
07965-SB	LOCATOR Implant Guided Drill, 2.4mmD x 16mmL



**SURGICAL INSTRUMENTS AND ACCESSORIES (CONTINUED)**


PART#	DESCRIPTION
07966-SB	LOCATOR Implant Guided Drill, 2.4mmD x 20mmL
07967-SB	LOCATOR Implant Guided Drill, 2.4mmD x 24mmL
07968-SB	LOCATOR Implant Guided Drill, 2.6mmD x 16mmL
07969-SB	LOCATOR Implant Guided Drill, 2.6mmD x 20mmL
07970-SB	LOCATOR Implant Guided Drill, 2.6mmD x 24mmL
07971-SB	LOCATOR Implant Guided Drill, 2.9mmD x 16mmL
07972-SB	LOCATOR Implant Guided Drill, 2.9mmD x 20mmL
07973-SB	LOCATOR Implant Guided Drill, 2.9mmD x 24mmL
07974-SB	LOCATOR Implant Guided Drill, 3.2mmD x 16mmL
07975-SB	LOCATOR Implant Guided Drill, 3.2mmD x 20mmL
07976-SB	LOCATOR Implant Guided Drill, 3.2mmD x 24mmL
07977-SB	LOCATOR Implant Guided Drill, 3.6mmD x 16mmL
07978-SB	LOCATOR Implant Guided Drill, 3.6mmD x 20mmL
07979-SB	LOCATOR Implant Guided Drill, 3.6mmD x 24mmL
07980-SB	LOCATOR Implant Guided Drill, 4.1mmD x 16mmL
07981-SB	LOCATOR Implant Guided Drill, 4.1mmD x 20mmL
07982-SB	LOCATOR Implant Guided Drill, 4.1mmD x 24mmL
07983-SB	LOCATOR Implant Guided Drill, 4.6mmD x 16mmL
07984-SB	LOCATOR Implant Guided Drill, 4.6mmD x 20mmL
07985-SB	LOCATOR Implant Guided Drill, 4.6mmD x 24mmL
07986-SB	LOCATOR Implant Guided Cortical Drill 3.3mm
07987-SB	LOCATOR Implant Guided Cortical Drill 3.7mm
07988-SB	LOCATOR Implant Guided Cortical Drill 4.3mm
07989-SB	LOCATOR Implant Guided Cortical Drill 4.8mm
07990-SB	LOCATOR Implant Guided Manual Implant Driver
07992-SB	LOCATOR Implant Guided Implant Latch Driver
07949	LOCATOR Implant Guided Surgical Tray Only

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-TX
LOCATOR CHAIRSIDE
MATERIALS & EQUIPMENT

# Attachment Systems **LOCATOR**<sup>®</sup>


**The self-aligning attachment innovation that revolutionized the industry and continues to be the market-leading attachment system.**

Engineered to flex, pivot and absorb, it's designed to protect what matters most. With over 20 years of clinical success, it's the most trusted overdenture attachment system manufactured only by Zest.




# LOCATOR® ATTACHMENT SYSTEM

LOCATOR®  
Scan Bodies  
\*not included with the abutment



(Processing Pack - 2 Pack)  
Standard Range or  
Extended Range

## STANDARD INSERT



Blue Low  
Retention



Pink Medium  
Retention



Clear High  
Retention

## EXTENDED RANGE INSERT



Gray Zero  
Retention



Red Low  
Retention



Orange Medium  
Retention



Green High  
Retention

## LOCATOR® PROCESSING COMPONENTS

PART#	DESCRIPTION
09530	Angle Measurement Guide
08517	Parallel Post - 4 Pack
08505	Impression Coping - 4 Pack
08505-20	Impression Coping - 20 Pack
08530	Abutment Analog (4mm Diameter) - 4 Pack
08530-20	Abutment Analog (4mm Diameter) - 20 Pack
08516	Abutment Analog (5mm Diameter) - 4 Pack
08516-20	Abutment Analog (5mm Diameter) - 20 Pack
08514	Block-Out Spacer - 20 Pack
08569	Processing Spacer - 4 Pack
08569-20	Processing Spacer - 20 Pack

LOCATOR Root

LOCATOR  
Implants

LOCATOR  
Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL  
HYGIENE

LOCATOR R-TX

LOCATOR  
CHAIRSIDE

MATERIALS &  
EQUIPMENT

## LOCATOR® PROCESSING COMPONENTS (CONTINUED)

PART#	DESCRIPTION
08519-2	LOCATOR Insert Processing Package - 2 Pack (Includes: Denture Cap with Processing Inserts, Block-Out Spacer, Standard Range Insert: Clear, Pink, Blue)
08519	LOCATOR Insert Processing Package - 1 Pack (Includes: Denture Cap with Processing Inserts, Block-Out Spacer, Standard Range Insert: Clear, Pink, Blue)
08519-10	LOCATOR Insert Processing Package - 10 Pack (Includes: Denture Cap with Processing Insert, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue)
08540	Extended Range Insert Processing Package - 1 Pack (Includes: Denture Cap with Processing Insert, Block-Out Spacer, Extended Range Inserts: Green, Orange, Red)
08540-2	Extended Range Insert Processing Package - 2 Pack (Includes: Denture Cap with Processing Insert, Block-Out Spacer, Extended Range Inserts: Green, Orange, Red)
08540-10	Extended Range Insert Processing Package - 10 Pack (Includes: Denture Cap with Processing Insert, Block-Out Spacer, Extended Range Inserts: Green, Orange, Red)
08550-2	LOCATOR Insert Processing Package SS (Stainless Steel) - 2 Pack (Includes: Stainless Steel Denture Cap with Processing Insert, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue)
08550-10	LOCATOR Insert Processing Package SS (Stainless Steel) - 10 Pack (Includes: Stainless Steel Denture Cap with Processing Insert, Block-Out Spacer, Standard Range Inserts: Clear, Pink, Blue)
08515	Black: Processing Replacement Insert - 4 Pack
08515-20	Black: Processing Replacement Insert - 20 Pack
08510	Replacement Denture Cap Assembly - 4 Pack
08510-10	Replacement Denture Cap Assembly - 10 Pack
08524	Clear: High Retention Replacement Insert - 4 Pack
08524-20	Clear: High Retention Replacement Insert - 20 Pack
08527	Pink: Medium Retention Replacement Insert - 4 Pack
08527-20	Pink: Medium Retention Replacement Insert - 20 Pack
08529	Blue: Low Retention Replacement Insert - 4 Pack
08529-20	Blue: Low Retention Replacement Insert - 20 Pack

## LOCATOR® PROCESSING COMPONENTS (CONTINUED)

PART#	DESCRIPTION
08547	Green: High Retention Extended Range Replacement Insert - 4 Pack
08547-20	Green: High Retention Extended Range Replacement Insert - 20 Pack
08915	Orange: Medium Retention Extended Range Replacement Insert - 4 Pack
08915-20	Orange: Medium Retention Extended Range Replacement Insert - 20 Pack
08548	Red: Low Retention Extended Range Replacement Insert - 4 Pack
08548-20	Red: Low Retention Extended Range Replacement Insert - 20 Pack
08558	Gray: Zero Retention Extended Range Replacement Insert - 4 Pack
08558-20	Gray: Zero Retention Extended Range Replacement Insert - 20 Pack
09260-02	LOCATOR Scan Body - 2 Pack
09260-04	LOCATOR Scan Body - 4 Pack
09260-10	LOCATOR Scan Body - 10 Pack

**NOTE: For Laboratory Processed Reline, Order Analogs and Processing Replacement Insert**

## LOCATOR® TOOLS

PART#	DESCRIPTION
08340	Enhanced LOCATOR Core Tool
08536	Manual LOCATOR & LOCATOR FIXED Driver
08393	LOCATOR Core Tool
08390	LOCATOR Abutment Driver (Gold End)
08317	LOCATOR Abutment Driver - Pitt Easy
08394	LOCATOR Abutment Holder Sleeve - 4 Pack
08394-20	LOCATOR Abutment Holder Sleeve - 20 Pack
08397	LOCATOR Core Tool, Insert Removal (Tip End)
09020	30Ncm Torque Wrench, Thumb Knob and 15mm Insert Driver
04391	20Ncm Torque Wrench, Thumb Knob and 15mm Insert Driver
08260	35Ncm Torque Wrench Insert Driver (15mm) for Straumann

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## LOCATOR® TOOLS (CONTINUED)

PART#	DESCRIPTION
08280	35Ncm Torque Wrench Insert Driver (21mm) for Straumann
08913	Torque Wrench Insert Driver (23mm), Latch Lock Connection
08926	LOCATOR Square Torque Wrench Insert Driver (15mm)
08927	LOCATOR Square Torque Wrench Insert Driver (21mm)
04936	.050 Hex Insert Driver

## LOCATOR® MULTI-UNIT ABUTMENTS (FREE STANDING APPLICATION)

PART#	DESCRIPTION
<b>NOTE: For current LOCATOR Multi-Unit Abutment Connections, please call 800.262.2310</b>	
08909-2	LOCATOR Abutment with Multi-Unit Titanium Collar - 2 Pack
08909-10	LOCATOR Abutment with Multi-Unit Titanium Collar - 10 Pack

## LOCATOR® PATIENT EDUCATION MODELS

PART#	DESCRIPTION
09015	LOCATOR Patient Education Model - 2 Implant Model

## LOCATOR® STARTER AND COMPONENT KITS

PART#	DESCRIPTION
09606-M	LOCATOR REMOVABLE Processing Component Storage Box (Box Only)
09605-P	LOCATOR REMOVABLE Processing Component Storage Organizer and Starter Kit Includes: Storage Box (1) Angle Measure Guide (1) LOCATOR Parallel Post 4 Pack (1) LOCATOR Scan Body 2 Pack (1) LOCATOR Impression Coping 4 Pack (1) LOCATOR Female Analog 4.0mm 4 Pack (1) LOCATOR Processing Spacer 4 Pack (1) LOCATOR Replacement Denture Cap Male Assembly 10 Pack (1) LOCATOR Extra Light Retention Replacement Male: Blue 4 Pack (1) LOCATOR Light Retention Replacement Male: Pink 4 Pack (1) LOCATOR Regular Retention Replacement Male Clear 4 Pack (1) LOCATOR Extended Range Extra Light Retention Replacement Male: Red 4 Pack (1) LOCATOR Extended Range Light Retention Replacement Male: Orange 4 Pack (1) LOCATOR Extended Range Regular Retention Replacement Male: Green 4 Pack (1) LOCATOR Zero Retention Male Extended Range Gray 4 Pack (1) LOCATOR Block Out Spacer FP 20 pack (1) LOCATOR Core Tool (1) Plastic Tweezers (1)



LOCATOR REMOVABLE Processing Component  
Storage Organizer Box - Box Only



LOCATOR REMOVABLE Processing Component  
Storage Organizer and Starter Kit



## LOCATOR® ABUTMENTS ARE IMPLANT SPECIFIC

TO ORDER, YOU WILL NEED THE FOLLOWING INFORMATION:

1. Name of the implant system
2. Implant diameter
3. Connection type and platform diameter
4. Tissue measurement (at the highest point)

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

### AB DENTAL

#### INTERNAL HEX CONNECTION

(Compatible with I2)

3.5mm, Narrow Platform, Blue	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm, Wide Platform, Black	08671	08672	08673	08674	08675	08627	08628

### ACE SURGICAL

#### INTERNAL HEX CONNECTION

(Compatible with Infinity)

3.5mm	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm	08671	08672	08673	08674	08675	08627	08628

#### TRI-CHANNEL

(Compatible with TRI-CAM)

3.5mm	08756	08757	08758	08759	08760	08734	08735
4.3mm	08761	08762	08763	08764	08765	08753	
5.0mm		08766	08767	08768	08769	08770	

#### TISSUE LEVEL

(Compatible with Infinity Octagon)

4.8mm, Regular Platform		08621 (1.25)	08622 (2.25)	08623 (3.25)	08624 (4.25)		
-------------------------	--	--------------	--------------	--------------	--------------	--	--

#### BONE LEVEL

(Compatible with Infinity)

2.7mm, Narrow Platform		02131	02132	02133	02134	02135	02136
3.3mm, Regular Platform		02137	02138	02139	02140	02141	02142

### ADIN IMPLANTS

#### INTERNAL HEX CONNECTION

(Compatible with Touareg, Swell)

3.5mm	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)		
-------	-------	-------	-------------	-------------	-------------	--	--

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## ALPHA BIO TECH

### INTERNAL HEX CONNECTION

(Compatible with Internal Hex)

3.75mm	01192(0.3)	01193	01194	01195	01196	01197	01198
--------	------------	-------	-------	-------	-------	-------	-------

## AMERICAN DENTAL IMPLANT CORP.

### INTERNAL HEX CONNECTION

(Compatible with Screw Vent)

3.5mm, Small Platform	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm, Wide Platform	08671	08672	08673	08674	08675	08627	08628
5.7mm, XL Platform		08180	08181	08182	08183	08184	08185

## ANTHOGRYR

### CONICAL CONNECTION

3.4mm Short	01697(1.5)	01698 (2.5)	01699 (3.5)	01016 (4.5)	01017 (5.5)		
-------------	------------	-------------	-------------	-------------	-------------	--	--

## BICON

### INTERNAL HEX CONNECTION

(Compatible with HA, Integra, Max)

2.0mm Post		08746	08747	08748	08749	08750	
2.5mm Post		02278	02279	02280	02281	02282	
3.0mm Post		08930	08931	08932	08933	08934	

## BIOHORIZONS

### INTERNAL HEX CONNECTION

(Compatible with Single Stage, Tapered, Tapered Plus)

3.5mm, Yellow	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm, Green	08671	08672	08673	08674	08675	08627	08628
5.0mm, Blue		08190	08191	08192	08193	08194	08195

(Compatible with Laser-Lok)

3.0mm, Gray		01927	01928	01929	01930	01931	01932
-------------	--	-------	-------	-------	-------	-------	-------

### EXTERNAL HEX CONNECTION

(Compatible with External Hex)

3.5mm, Yellow		08531 (1.25)	08532	08533	08534	08535	
4.0mm, Green		08541	08542	08543	08544	08545	
5.7mm, Blue		08551 (1.5)	08552 (2.5)	08553 (3.5)	08554 (4.5)	08555 (5.5)	



IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## BIOHORIZONS (continued)

### TAPERED PRO CONICAL CONNECTION



(Compatible with CONELOG)

3.3 mm Narrow Platform		02748	02749	02750	02751	02752	02753
3.8 mm Regular Platform		02754	02755	02756	02757	02758	02759

## BIO-LOK

### EXTERNAL HEX CONNECTION

(Compatible with Micro-Lok, Silhouette)

3.45mm, Green		08561	08562	08563	08564	08565	
4.0mm, Yellow	08681 (0.73)		08682	08683	08684	08685	08696
5.0mm, Red		08581 (1.4)	08582	08583	08584	08585	

## BIOMET 3i

SEE ZIMVIE SECTION

## BLUE SKY BIO

### INTERNAL HEX CONNECTION

(Compatible BIO-Internal Hex)

3.5mm	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm	08671	08672	08673	08674	08675	08627	08628

### TRI-CHANNEL

(Compatible with TRILOBE)

3.5mm, Narrow Platform, Purple	08756	08757	08758	08759	08760	08734	08735
4.3mm, Regular Platform, Gold	08761	08762	08763	08764	08765	08753	08754
5.0mm/6.0mm, Wide Platform, Blue		08766	08767	08768	08769	08770	

### TISSUE LEVEL

(Compatible with One Stage)

4.8mm, Regular Platform	08950 (0.73)		08951	08952	08953	08954	08955
-------------------------	--------------	--	-------	-------	-------	-------	-------

### BONE LEVEL

(Compatible with BIO Quattro)

2.8mm, Narrow Platform		02131	02132	02133	02134	02135	02136
3.3mm, Regular Platform		02137	02138	02139	02140	02141	02142

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## BTI-BIOTECHNOLOGY INSTITUTE SL

### EXTERNAL HEX CONNECTION

(Compatible with Externa, Externa Tiny)

3.5mm, Narrow Platform		02979	02980	02981	02982	02983	02984
4.1mm, Standard Platform	08681 (0.73)		08682	08683	08684	08685	08696

### INTERNAL HEX CONNECTION

(Compatible with Interna)

4.1mm, Standard Platform, Blue/Yellow		08291	08292	08293	08294	08295	08296
---------------------------------------	--	-------	-------	-------	-------	-------	-------

## DENTIUM

### CONICAL CONNECTION

(Compatible with Implantium and Superline)

4.0mm	08656(0.5)		08657	08658	08659	08660	
-------	------------	--	-------	-------	-------	-------	--

## DENTSPLY

### CONICAL CONNECTION

(Compatible with Astra Tech EV)

3.6mm		02790	02791	02792	02793	02794	
4.2mm		02795	02796	02797	02798	02799	
4.8mm		02800	02801	02802	02803	02804	
2.9mm		08115	08116	08117	08118	08119	
3.9mm	08686 (0.5)		08687	08688	08689	08690	

### CONICAL CONNECTION

(Compatible with Osseospeed TX)

3.5/4.0mm	08115 (.81)		08116	08117	08118	08119	
4.5/5/0mm	08656(0.5)		08660	08661	08662	08663	

### INTERNAL HEX CONNECTION

(Compatible with FRIALIT, XIVE)

3.4mm		08805	08806	08807	08808	08809	
3.8mm		08810	08811	08812	08813	08814	
4.5mm		08815	08816	08817	08818	08819	
5.5mm		08139 (1.3)	08140	08141	08142	08143	

<b>IMPLANT COMPANY</b>	<b>LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)</b>						
<b>DESCRIPTION</b>	<b>CUFF HEIGHT</b>						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## DIO USA

### CONICAL CONNECTION

**(Compatible with SM Submerged)**

3.1mm, Narrow Platform		01351	01352	01353	01354	01355	01356
3.5mm, Regular/Wide Platform		01361	01362	01363	01364	01365	01366

**(Compatible with UF(II))**

2.3mm, Narrow Platform		02961	02962	02963	02964	02965	02966
3.35mm, Regular Platform		02967	02968	02969	02970	02971	02972

### GLIDEWELL/PRISMATIK DENTALCRAFT

#### INTERNAL HEX CONNECTION

**(Compatible with Inclusive)**

3.0mm, Aqua		02212	02213	02214	02215	02216	02217
3.5mm, Green	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm, Purple	08671	08672	08673	08674	08675	08627	08628

### GLIDEWELL/PRISMATIK DENTALCRAFT

#### INTERNAL HEX CONNECTION

**(Compatible with HAHN)**

3.0mm		02212	02213	02214	02215	02216	02217
3.5mm		02069	02070	02071	02072	02073	02074
5.0mm		02075	02076	02077	02078	02079	02080

### HIOSSEN/OSSTEM

#### INTERNAL HEX CONNECTION

**(Compatible with EK)**

3.5 mm, Single		02392-SB	02393-SB	02394-SB	02395-SB	02396-SB	02397-SB
----------------	--	----------	----------	----------	----------	----------	----------

#### INTERNAL HEX CONNECTION

**(Compatible with ET/TS)**

2.8mm, Mini Platform		01561	01562	01563	01564	01565	01566
3.35mm, Standard Platform		01551	01552	01553	01554	01555	01556

#### INTERNAL HEX CONNECTION

**(Compatible with SS)**

4.8mm, Regular Platform	08950 (0.73)		08951	08952	08953	08954	08955
-------------------------	--------------	--	-------	-------	-------	-------	-------

LOCATOR Root  
LOCATOR Implants  
LOCATOR Abutments  
LOCATOR Bar  
LOCATOR FIXED  
REGAL HYGIENE  
LOCATOR R-TX  
LOCATOR CHAIRSIDE  
MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## HIOSEN/OSSTEM (continued)

### EXTERNAL HEX CONNECTION

#### (Compatible with US)

3.5mm, Mini Platform	08105 (0.73)		08106	08107	08108	08109	
4.1mm, Regular Platform	08681 (0.73)		08682	08683	08684	08685	08696
5.1mm, Wide Platform		08825 (1.1)	08826	08827	08828	08829	

## IBS OF AMERICA

### CONICAL CONNECTION

#### (Compatible with Anker, N.R FIX, MAGIC FC)

3.0mm		01551	01552	01553	01554	01555	01556
-------	--	-------	-------	-------	-------	-------	-------

## IMPLANT DIRECT

### INTERNAL HEX CONNECTION

#### (Compatible with Legacy)

3.0mm, Blue	02308 (0.75)		02309	02310	02311	02312	
3.5mm, Green	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm, Purple	08671	08672	08673	08674	08675	08627	08628
5.7mm, Yellow		01471	01472	01473	01474	01475	01476

### INTERNAL HEX CONNECTION

#### (Compatible with Innova-Endopore)

4.1mm	08145	08146	08147	08148	08149	08150	
5.0mm		08164	08165	08166	08167	08168	08169

#### (Compatible with Pitt-Easy)

3.3mm		08311	08312	08313	08314	08315	08316
3.75mm		08321	08322	08323	08324	08325	08326
4.1mm		08331	08332	08333	08334	08335	08336
4.9mm		08341	08342	08343	08344	08345	08346

#### (Compatible with SID)

3.7mm			01250	01251	01252	01253	01254
-------	--	--	-------	-------	-------	-------	-------

### CONICAL CONNECTION

#### (Compatible with Swish-Active, Interactive, Simply Iconic)

3.0mm, Narrow Platform, Purple		02069	02070	02071	02072	02073	02074
-----------------------------------	--	-------	-------	-------	-------	-------	-------

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

### IMPLANT DIRECT (continued)

#### CONICAL CONNECTION

3.4mm, Regular Platform, Yellow		02075	02076	02077	02078	02079	02080
---------------------------------	--	-------	-------	-------	-------	-------	-------

#### TISSUE LEVEL CONNECTION

**(Compatible with Swish-Plus, Swish-Tapered)**

3.7mm, Blue	02308 (0.75)		02309	02310	02311	02312	02313
4.8mm, Yellow		02106	02107	02108	02109	02110	02111

#### TRI-CHANNEL CONNECTION

**(Compatible with Replant, Reactive, Replus)**

3.5mm, Purple	08756	08757	08758	08759	08760	08734	08735
4.3mm, Yellow	08761	08762	08763	08764	08765	08753	08754
5.0mm, Blue		08766	08767	08768	08769	08770	

#### EXTERNAL HEX CONNECTION

**(Compatible with Screwplant, Screwplus)**

3.7mm, Green		01371	01372	01373	01374	01375	01376
4.7mm, Purple		01381 (1.5)	01382	01383	01384	01385	01386
5.7mm, Yellow		01391	01392	01393	01394	01395	01396

**(Compatible with SwissPlant)**

3.7mm	01821 (0.4)	01822	01823	01824	01825	01826	
-------	-------------	-------	-------	-------	-------	-------	--

#### EXTERNAL HEX CONNECTION

**(Compatible with Sybron-Endopore)**

3.3mm		01631	01632	01633	01634	01635	01636
4.1mm		01637	01638	01639	01640	01641	01642

**(Compatible with Innova- Sybron)**

4.1mm	08681 (0.73)		08682	08683	08684	08685	08686
-------	--------------	--	-------	-------	-------	-------	-------

**(Compatible with Sybron Ex Hex Plus)**

3.5mm	08105 (0.73)		08106	08107	08108	08109	
4.1mm	08681 (0.73)		08682	08683	08684	08685	08696
5.0mm		08825	08826	08827	08828	08829	

**(Compatible with Innova- Endopore)**

3.5mm		08851	08852	08853	08854	08855	
4.8mm		08950	08951	08952	08953	08954	08955
5.0mm		08691 (1.35)	08692	08693	08694	08695	

LOCATOR Root  
LOCATOR Implants  
LOCATOR Abutments  
LOCATOR Bar  
LOCATOR FIXED  
REGAL HYGIENE  
LOCATOR R-TX  
LOCATOR CHAIRSIDE  
MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## IMPLANT LOGISTICS



**(Compatible with Implant One 300 Series)**

2.75mm		02180	02181	02182	02183	02184	02185
--------	--	-------	-------	-------	-------	-------	-------

**(Compatible with Implant One 400 Series)**

3.25mm		02187	02188	02189	02190	02191	02192
--------	--	-------	-------	-------	-------	-------	-------

### 3M

#### INTERNAL HEX CONNECTION

**(Compatible with Endure)**

4.3mm	08994 (0.63)		08995	08996	08997	08998	
-------	-----------------	--	-------	-------	-------	-------	--

#### INTERNAL HEX CONNECTION

**(Compatible with Mini, Universal)**

3.3mm	08856 (0.91)		08857	08858	08859	08860	
3.4mm/4.0mm	08681 (0.73)		08682	08683	08684	08685	

### IMZ

#### NON-HEX CONNECTION (FLAT TOP)

3.3mm	08020	08021	08022	08023	08024	08025	
4.0mm	08030 (0.35)	08031	08032	08033	08034	08035	

#### EXTERNAL HEX CONNECTION

3.3mm/4.0mm	08681 (0.73)		08682	08683	08684	08685	
-------------	-----------------	--	-------	-------	-------	-------	--

### INTRA-LOCK

#### EXTERNAL HEX CONNECTION

**(Compatible with External Hex)**

4.1mm	08681 (0.73)		08682	08683	08684	08685	08696
5.0mm		08727 (1.42)	08728	08729	08730	08731	

**(Compatible with External Hex Conic)**

4.3mm		08781	08782	08783	08784	08785	
5.0mm		08727 (1.42)	08728	08729	08730	08731	

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## IQ IMPLANTS

### INTERNAL HEX CONNECTION

(Compatible with Smart and Master)

3.75mm	01192 (0.3)	01193	01194	01195	01196	01197	01198
--------	----------------	-------	-------	-------	-------	-------	-------

## KAT IMPLANTS

### MORSE TAPER CONNECTION

(Compatible with KAT Implants)

4.0mm			01714 (2.5)	01715 (3.5)	01716 (4.5)		
-------	--	--	----------------	----------------	----------------	--	--

## KEYSTONE DENTAL

### TILOBE CONNECTION

(Compatible with Genesis, Primaconnex, Max TI)

3.5mm		08073	08074	08075	08076	08077	08097
4.5mm		08079	08080	08081	08082	08083	08098
5.0mm		08085	08086	08087	08088	08089	08099

## KEYSTONE DENTAL

### INTERNAL HEX CONNECTION

(Compatible with Renova)

3.75mm, SDI	08090	08091	08092	08093	08094	08095	08096
4.75mm, RDI	08125 (0.63)		08126	08127	08128	08129	08130

(Compatible with Sustain)

4.0mm	08651	08652	08653	08654	08655		
-------	-------	-------	-------	-------	-------	--	--

### MORSE TAPER CONNECTION

(Compatible with Stage-1, OCT-XT)

4.8mm	08950		08951	08952	08953	08954	08955
-------	-------	--	-------	-------	-------	-------	-------

### EXTERNAL HEX CONNECTION

(Compatible with Restore, Sustain)

3.3mm, SD, Green		08196	08197	08198	08199	08200	
4.1mm, RD, Blue	08681 (0.73)		08682	08683	08684	08685	08696
5.0mm, WD, Purple	08727 (1.42)		08728	08729	08730	08731	

### EXTERNAL HEX CONNECTION

(Compatible with TRI-NEX)

3.5mm, Purple	08756	08757	08758	08759	08760	08734	08735
4.3mm, Yellow	08761	08762	08763	08764	08765	08753	08754
5.0mm, Blue		08766	08767	08768	08769	08770	

LOCATOR Root  
LOCATOR Implants  
LOCATOR Abutments  
LOCATOR Bar  
LOCATOR FIXED  
REGAL HYGIENE  
LOCATOR R-TX  
LOCATOR CHAIRSIDE  
MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

KYOCERA							
BUTT JOINT CONNECTION							
4.2mm	08716 (0.35)	08717	08718	08719	08720	08721	

LASAK							
BUTT JOINT CONNECTION							
3.7mm	01207	01208	01209	01210	01211	01212	

MEGAGEN							
INTERNAL HEX CONNECTION							
<b>(Compatible with AnyRidge)</b>							
3.5mm		02100	02101	02102	02103	02104	02105

INTERNAL OCTA CONNECTION							
<b>(Compatible with Blue Diamond)</b>							
2.8 mm		02131	02132	02133	02134	02135	02136
3.3mm		02406	02407	02408	02409	02410	02411

MIS							
INTERNAL HEX CONNECTION							
<b>(Compatible with Seven, M4)</b>							
3.3mm, Narrow Platform, Yellow	02062	02063	02064	02065	02066	02067	02068
3.75mm, Standard Platform, Purple	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm, Wide Platform, Green	08671	08672	08673	08674	08675	08627	08628

EXTERNAL HEX CONNECTION							
<b>(Compatible with Biocom, Lance)</b>							
5.0mm, Wide Platform	08820 (0.6)		08821	08822	08823	08824	
3.75/4.2mm Standard Platform	08681(.73)		08682	08683	08684	08685	08696

CONICAL CONNECTION							
<b>(Compatible with C1, V3)</b>							
2.75mm, Narrow Platform, Yellow/Blue		02143	02144	02145	02146	02147	02148
3.15mm, Standard Platform, Purple		02149	02150	02151	02152	02153	02154



IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## NEODENT

### CONE MORSE/GRAND MORSE CONNECTION

**NEW!**

#### (Compatible with GM)

2.9 mm Narrow		02173	02174	02175	02176	02177	02178
3.0mm		02212	02213	02214	02215	02216	02217

#### (Compatible with CM)

2.0mm			02168	02169	021610	021611	021612
-------	--	--	-------	-------	--------	--------	--------

## NEOSS

### INTERNAL HEX CONNECTION

#### (Compatible with ProActive Straight, ProActive Tapered)

4.0mm		08451	08452	08453	08454	08455	08456
-------	--	-------	-------	-------	-------	-------	-------

## NOBEL BIOCARE

### INTERNAL HEX CONNECTION

#### (Compatible with Steri-Oss)

3.25mm	08611	08612	08613	08614	08615		
4.0mm	08651	08652	08653	08654	08655		

### EXTERNAL HEX CONNECTION

#### (Compatible with NobelReplace, Steri-Oss)

3.5mm, Narrow Platform		08676	08677	08678	08679	08680	
4.0mm, Standard Platform		08781	08782	08783	08784	08785	
5.0mm, Wide Platform		08835 (1.37)	08836	08837	08838	08839	

#### (Compatible with Hex-Lock)

3.25mm		08646 (1.5)	08647 (2.5)	08648 (3.5)	08649 (4.5)	08650 (5.5)	
4.0mm	08681 (0.73)		08682	08683	08684	08685	08696
5.0mm		08835 (1.37)	08836	08837	08838	08839	

#### (Compatible with Branemark, NobelSpeedy Groovy)

3.5mm, Narrow Platform, Grey	08105		08106	08107	08108	08109	
4.1mm, Regular Platform, Grey	08681 (0.73)		08682	08683	08684	08685	
5.1mm, Wide Platform, Grey		08825 (1.1)	08826	08827	08828	08829	

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## NOBEL BIOCARE (continued)

### CONICAL CONNECTION

(Compatible with Nobel Active)

3.5mm, Narrow Platform, Purple		02069	02070	02071	02072	02073	02074
4.3mm, Regular Platform, Yellow		02075	02076	02077	02078	02079	02080

### TRI-CHANNEL CONNECTION

(Compatible with Nobel Replace)

3.5mm, Narrow Platform, Purple	08756	08757	08758	08759	08760	08734	08735
4.3mm, Regular Platform, Yellow	08761 (0.25)	08762	08763	08764	08765	08753	08754
5.0mm, Wide Platform, Blue		08766	08767	08768	08769	08770	

## OCO BIOMEDICAL

### INTERNAL HEX CONNECTION

(Compatible with Engage)

3.0mm	01707	01708	01709	01710	01711	01712	01713
3.5mm	08661	08662	08663	08664	08665	08625	08626
4.5mm	08671	08672	08673	08674	08675	08627	
6.0mm		01471	01472 (2.5)	01473 (3.5)	01474 (4.5)	01475 (5.5)	01476 (6.5)

## OSSOTANIUM

### CONICAL CONNECTION

(Compatible with Basic, Osso-tight, BL)

3.0mm, Narrow Platform		02069	02070	02071	02072	02073	02074
3.4mm, Regular Platform		02075	02076	02077	02078	02079	02080

## OSTEOREADY

### INTERNAL HEX CONNECTION

(Compatible with Seven, M4)

3.5mm	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)	08625 (5.5)	08626 (6.5)
4.5mm	08671	08672	08673	08674	08675	08627	08628

## PALTOP

### INTERNAL HEX CONNECTION

(Compatible with Advanced)

2.9mm, Narrow Platform		02296	02297	02298	02299	02300	02301
3.65mm, Standard Platform		02284	02285	02286	02287	02288	02289

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## PALTOP (continued)

### CONICAL CONNECTION

(Compatible with Active)

2.9mm		02290	02291	02292	02293	02294	02295
-------	--	-------	-------	-------	-------	-------	-------

### INTERNAL HEX CONNECTION

(Compatible with TSA)

3.7mm		01755	01756	01757	01758	01759	01760
4.7mm		01761	01762	01763	01764	01765	01766
6.0mm			01767	01768	01769	01770	01771

### EXTERNAL HEX CONNECTION

(Compatible with TSH, BNT)

3.3mm		01916	01917	01918	01919	01920	01921
4.0mm		01942	01943	01944	01945	01946	01947
5.0mm			01922	01923	01924	01925	01926

## RITTER IMPLANTS

### INTERNAL HEX CONNECTION

(Compatible with SB/LA Spiral, Quadro Spiral)

3.0mm, Narrow Platform	02308 (0.75)		02309	02310	02311	02312	02313
------------------------	-----------------	--	-------	-------	-------	-------	-------

## SARGON

### External Hex Connection

3.8mm		08722		08723	08724	08725	08726
-------	--	-------	--	-------	-------	-------	-------

## SHINHUNG

### CONICAL CONNECTION

(Compatible with LUNA)

2.8mm, Narrow Platform		01561	01562	01563	01564	01565	01566
3.35mm, Regular Platform		01551	01552	01553	01554	01555	01556

## SOUTHERN IMPLANTS

### INTERNAL OCTAGON CONNECTION

(Compatible with ITC, ITT)

4.8mm	08950 (0.73)		08951	08952	08953	08954	08955
3.5mm	08756	08757	08758	08759	08760	08734	08735
4.3mm	08761	08762	08763	08764	08765	08753	08754
5.0mm		08766	08767	08768	08769	08770	

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-TX
LOCATOR CHAIRSIDE
MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## STRAUMANN

### TISSUE LEVEL CONNECTION

**(Compatible with Roxolid)**

3.3mm, NNC		02056	02057	02058	02059	02060	02061
4.1mm, RN		08742	08743	08744	08745	08751	08752
4.8mm, WN		08870	08871	08872	08873	08874	

**(Compatible with TLX)**

NT TLX		02350-SB	02351-SB	02352-SB	02353-SB	02354-SB	02355-SB
RT TLX		02356-SB	02357-SB	02358-SB	02359-SB	02360-SB	02361-SB
WT TLX		02362-SB	02363-SB	02364-SB	02365-SB	02366-SB	02367-SB

### BONE LEVEL CONNECTION

**(Compatible with Roxolid)**

3.3mm, NC			08422	08423	08424	08425	08426
4.1mm, RC		08170	08171	08172	08173	08174	08175

**(Compatible with BLX)**

3.5mm, RB/WB		02326	02327	02328	02329	02330	02331
--------------	--	-------	-------	-------	-------	-------	-------

## TATUM

### INTERNAL HEX CONNECTION

**(Compatible with Tapered Internal)**

3.5mm	08100	08101	08102	08103	08104		
4.0mm	08110 (0.25)	08111	08112	08113	08114		
4.5mm	08120 (0.5)	08121	08122	08123	08124		

### CONICAL CONNECTION

**(Compatible with Integrity)**

3.25mm		02314	02315	02316	02317	02318	02319
4.0mm		02320	02321	02322	02323	02324	02325
5.0mm		02155	02156	02157	02158	02159	02160

<b>IMPLANT COMPANY</b>	<b>LOCATOR® ATTACHMENT SYSTEM (Includes Abutment)</b>						
<b>DESCRIPTION</b>	<b>CUFF HEIGHT</b>						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## ZimVie/ BIOMET 3i

### INTERNAL CONNECTION

**(Compatible with Certain, PREVAIL, and Zimvie T3 PRO)**

3.4mm		08431	08432	08433	08434	08435	08436
4.1mm		08630	08631	08632	08633	08634	08635

### INTERNAL HEX CONNECTION

**(Compatible with Screw Vent, Tapered Screw Vent, Micro/Core/Bio Vent)**

3.5mm, Green	08661	08662	08663 (2.5)	08664 (3.5)	08665 (4.5)		
4.5mm, Purple	08671	08672	08673	08674	08675		
5.7mm, Yellow		08180	08181	08182	08183	08184	08185

**(Compatible with Advent)**

4.5mm, Purple	08841 (0.71)		08842	08843	08844	08845	08846
---------------	-----------------	--	-------	-------	-------	-------	-------

### CONICAL CONNECTION

**(Compatible with Eztetic)**

2.9mm		02218	02219	02220	02221	02222	02223
-------	--	-------	-------	-------	-------	-------	-------

**(Compatible with Spline)**

3.25mm			08701	08702	08703	08704	08705
4.0mm		08711	08712	08713	08714	08715	
5.0mm			08736	08737	08738	08739	08740

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE


LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT



## LOCATOR® BAR ATTACHMENT



PART#	DESCRIPTION
<b>All LOCATOR Bars, Including Multi-Unit, Require 08028 LOCATOR Yellow Bar Male Processing Package</b>	
08013	Bar 2-56 Castable Threaded Insert - 10 Pack
08587-2	Bar, 2-56 Thread - 2 Pack
08587-10	Bar, 2-56 Thread - 10 Pack
08014	Bar 2.0mm Castable Threaded Insert - 10 Pack
08589-2	Bar, 2.0mm Thread - 2 Pack
08589-10	Bar, 2.0mm Thread - 10 Pack
08588-2	Laser Bar, Titanium - 2 Pack
08588-10	Laser Bar, Titanium - 10 Pack
08590-2	Laser Bar, Stainless Steel - 2 Pack
08590-10	Laser Bar, Stainless Steel - 10 Pack
08586-2	Cast-To Bar - 2 Pack
08586-10	Cast-To Bar - 10 Pack

LOCATOR Root

LOCATOR  
Implants

LOCATOR  
Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL  
HYGIENE

LOCATOR R-TX

LOCATOR  
CHAIRSIDE

MATERIALS &  
EQUIPMENT

## LOCATOR® FOR MULTI-UNIT ABUTMENT (BAR-SPLINTED APPLICATION)

**All LOCATOR Bars, Including Multi-Unit,  
Require 08028 LOCATOR Yellow Bar Insert Processing Package**

PART#	DESCRIPTION
08917-2	LOCATOR Abutment w/Castable Collar - 2 Pack
08917-10	LOCATOR Abutment w/Castable Collar - 10 Pack

### LOCATOR BAR INSERT PROCESSING PACKAGES

PART#	DESCRIPTION
08028-2	Bar Insert Processing Package - 2 Pack (Includes Denture Cap with Yellow Processing Insert, Standard Range Inserts: Clear, Pink, Blue and Block-Out Spacer)
08028-10	Bar Insert Processing Package - 10 Pack (Includes Denture Cap with Yellow Processing Insert, Standard Range Inserts: Clear, Pink, Blue and Block-Out Spacer)
08026	Yellow: Bar Processing Replacement Insert - 4 Pack
08026-20	Yellow: Bar Processing Replacement Insert - 20 Pack

### LOCATOR BAR ATTACHMENT TOOLS

PART#	DESCRIPTION
08393	LOCATOR Core Tool
09103	1.8mm Bar Drill, 2-56 Thread
08340	Enhanced LOCATOR Core Tool
08536	Manual LOCATOR & LOCATOR FIXED Driver
09105	2-56 Bar Tap, 2-56 Thread
09102	1.7mm Bar Drill, 2.0mm Thread
09104	2.0mm Bar Tap, 2.0mm Thread
08016	Bar Drill & Tap Holder Tool
09107	Paralleling Mandrel




# Attachment Systems

# LOCATOR

# FIXED<sup>®</sup>


**Introducing a NEW flexible and accessible FIXED full-arch solution.**

- LOCATOR Removable or FIXED... it's your choice.
- Utilize patients' existing LOCATOR Abutments to transition to a FIXED full-arch solution.
- Offer staged treatment planning for your patients, removable today... FIXED tomorrow.
- Give your patients the life they had with their permanent teeth with an affordable, permanent teeth solution
- A LOCATOR proven alternative to the discomfort and pain associated with conventional dentures that feels more like permanent teeth.




## LOCATOR FIXED® ATTACHMENT SYSTEM

### LOCATOR FIXED® Abutment



### Inserts & Housing

Four Unit FIXED Insert




Green

Anterior/Posterior  
FIXED Insert



Tan

Mid-Arch FIXED Insert



Blue

FIXED Housing Assembly



Seating and Removal Tool



LOCATOR® Systems  
Enhanced Core Tool



## LOCATOR FIXED STARTER AND COMPONENT KITS

PART#	DESCRIPTION
09601-P	<b>LOCATOR FIXED Bronze Starter Package</b> Includes: LOCATOR Systems Enhanced CORE Tool (1); LOCATOR FIXED Seating and Removal Tool (1); LOCATOR FIXED Insert Anterior/Posterior, Tan, 2 Pack (2); LOCATOR FIXED Insert Mid Arch, Blue, 2 Pack (1); LOCATOR FIXED Insert 4 Unit, Green, 2 Pack (2); LOCATOR FIXED Housing Assembly, 4 Pack (1); LOCATOR FIXED Processing Pack, 4 Unit, Green (12); CHAIRSIDE Trepphine Bur (1)
09602-P	<b>LOCATOR FIXED Silver Starter Package</b> Includes: LOCATOR Systems Enhanced CORE Tool (1); LOCATOR FIXED Seating and Removal Tool (1); LOCATOR FIXED Insert Anterior/Posterior, Tan, 2 Pack (2); LOCATOR FIXED Insert Mid Arch, Blue, 2 Pack (1); LOCATOR FIXED Insert 4 Unit, Green, 2 Pack (2); LOCATOR FIXED Housing Assembly, 4 Pack (1); LOCATOR FIXED Processing Pack, 4 Unit, Green (12); LOCATOR FIXED Processing Pack, Anterior/Posterior Tan (6); LOCATOR FIXED Processing Pack Mid Arch Blue (4)
09575-P	<b>LOCATOR FIXED Gold Starter Package</b> Includes: LOCATOR Systems Enhanced CORE Tool (1); LOCATOR FIXED Seating and Removal Tool (1); LOCATOR FIXED Insert Anterior/Posterior, Tan, 2 Pack (2); LOCATOR FIXED Insert Mid Arch, Blue, 2 Pack (1); LOCATOR FIXED Insert 4 Unit, Green, 2 Pack (2); LOCATOR FIXED Housing Assembly, 4 Pack (1); LOCATOR FIXED Processing Pack, 4 Unit, Green (16); LOCATOR FIXED Processing Pack, Anterior/Posterior Tan (9); LOCATOR FIXED Processing Pack Mid Arch Blue (6)
09264-P	<b>LOCATOR FIXED Lab Starter Package</b> Includes: LOCATOR FIXED Housing Assembly, 10 Pack (3); LOCATOR FIXED® Insert Anterior/Posterior, 10 Pack (1); LOCATOR FIXED® Insert Mid Arch, Blue, 10 Pack (1); LOCATOR FIXED® Insert 4 Unit, Green, 10 Pack (1)
09263-P	<b>LOCATOR FIXED Tooling Package</b> Includes: 1 Core Tool and 1 Seating and Removal Tool

## LOCATOR FIXED STARTER AND COMPONENT KITS (CONTINUED)



LOCATOR FIXED Processing Component Storage Organizer Box - Box Only



LOCATOR FIXED Processing Component Storage Organizer and Starter Kit

PART#	DESCRIPTION
09604-M	LOCATOR FIXED Processing Storage Kit (Box Only)
09603-P	<b>LOCATOR FIXED Processing Storage Kit and Starter Kit</b> Includes: 1 Storage Box (1)Angle Measure Guide (1)LOCATOR Parallel Post 4 Pack (1) LOCATOR Scan Body 2 Pack (2)LOCATOR Impression Coping 4 Pack (1) LOCATOR Female Analog 4.0mm 4 Pack (1)LOCATOR Processing Spacer 4 Pack (1) LOCATOR FIXED Housing 10 Pack (1)LOCATOR FIXED Insert Anterior/Posterior, Tan, 2 Pack (2) LOCATOR FIXED Insert Mid Arch, Blue, 2 Pack (2)LOCATOR FIXED Insert 4 Unit, Green, 2 Pack (2)LOCATOR Black Processing Inserts 4 Pack (1)LOCATOR 35NCM Torque Wrench Drives 15MM (1)LOCATOR Block Out Spacers 20 Pack (1) LOCATOR Systems Core Tool (1)

## LOCATOR FIXED INSERTS AND HOUSINGS

PART#	DESCRIPTION
09506-2-SB	LOCATOR FIXED Insert Anterior/Posterior, Tan, 2 Pack
09507-2-SB	LOCATOR FIXED Insert Mid Arch, Blue, 2 Pack
09501-2-SB	LOCATOR FIXED Insert 4 Unit, Green, 2 Pack
09506-10-SB	LOCATOR FIXED Insert Anterior/Posterior, 10 Pack
09507-10-SB	LOCATOR FIXED Insert Mid Arch, Blue, 10 Pack
09501-10-SB	LOCATOR FIXED Insert 4 Unit, Green, 10 Pack
09508-4-SB	LOCATOR FIXED Housing Assembly, 4 Pack
09508-10-SB	LOCATOR FIXED Housing Assembly, 10 Pack
09509-01-SB	LOCATOR FIXED Processing Pack, 4 Unit, Green
09509-02-SB	LOCATOR FIXED Processing Pack, Anterior/Posterior Tan
09509-03-SB	LOCATOR FIXED Processing Pack Mid Arch Blue

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## LOCATOR FIXED TOOLING & OTHER

PART#	DESCRIPTION
08340-M	LOCATOR Systems Enhanced CORE Tool
08536-SB	LOCATOR Manual Abutment Driver
09599-M	LOCATOR FIXED Seating and Removal Tool
09260-2-SB	LOCATOR Scan Body, 2-Pack
09260-4-SB	LOCATOR Scan Body, 4 Pack
09260-10-SB	LOCATOR Scan Body, 10 Pack
08301-SB	LOCATOR Systems Core Tool Replacement O-Rings (4 pk)
08348-SB	LOCATOR Systems Core Tool Collet
10157-01-SB	LOCATOR FIXED Seating Tip
10156-01-M	LOCATOR FIXED Removal Tip
10095-SB	LOCATOR FIXED Optional Removal Hook Tip
10159-01-SB	LOCATOR FIXED Removal Tip Wire Loop
10158-01-SB	LOCATOR FIXED Seating Tip Cushion
10154-01-SB	LOCATOR FIXED Tool 2.5mm Hex Wrench
10155-01-SB	LOCATOR FIXED Tool Tip Wrench

## LOCATOR FIXED PATIENT EDUCATION

PART#	DESCRIPTION
L8084	LOCATOR FIXED Patient Brochure (25 Pack)
09166-M	LOCATOR FIXED Patient Education Model
L8086	LOCATOR FIXED Patient Testimonial Digital Book
L8087	LOCATOR FIXED Digital Patient Marketing Materials USB
L8092	LOCATOR FIXED Patient Posters (3 posters)
L8091	LOCATOR FIXED Patient Education Flip Chart
09610-M	LOCATOR FIXED Patient Marketing Starter Kit



# An Innovative Oral Care System

# REGAL HYGIENE

## Go and Protect!

Zest Dental Solutions, proud owner of the renowned LOCATOR Removable Overdenture and LOCATOR FIXED Permanent Teeth Systems, is excited to Announce...Regal Hygiene!



### Regal Protect:

- Lasts six months to one year.
- Prevents plaque & bacteria from adhering to teeth, dentures & prostheses.
- Resists acids, sticking foods, and microbiota.
- Clear, tasteless, and odorless silicone polymer gel coating.

### Regal GO:

- Removes debris and biofilm for a refreshing smile.
- Easy and effective cleaning solution for FIXED full-arch patients.
- The advanced brush design fits into hard-to-reach areas around implants, orthodontic appliances, and under prostheses.
- Regal GO contains EDTA, which disrupts biofilm and natural sweeteners for a fresher clean-mouth sensation.

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## REGAL HYGIENE


PART#	DESCRIPTION
90-00239-M	Regal GO, 10pk (10x4pk, 4ml each)
95614-M	Regal Protect, 5pk
90-00242-M	Regal Protect & GO Bundle (5 vials +5 mini Regal Go pens +5 patient brochures)



# Attachment Systems **LOCATOR** **R-Tx<sup>®</sup>**

**Durability and confidence that lasts.**

Redefined, re-imagined and improved.  
A better, simpler, stronger removable  
overdenture attachment innovation.



LOCATOR Root

LOCATOR  
Implants

LOCATOR  
Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL  
HYGIENE

LOCATOR R-TX


LOCATOR  
CHAIRSIDE

MATERIALS &  
EQUIPMENT

## LOCATOR R-Tx® REMOVABLE ATTACHMENT SYSTEM



### Standard Range Retention Inserts



### Limited Range Retention Inserts



## LOCATOR R-Tx® PROCESSING COMPONENTS


PART#	DESCRIPTION
30001-01	Gray: Zero Retention Insert - 4 Pack
30002-01	Blue: Low Retention Insert - 4 Pack
30003-01	Pink: Medium Retention Insert - 4 Pack
30004-01	Clear: High Retention Insert - 4 Pack
30027-01	Aqua: Low Limited Range Retention Insert - 4 Pack
30028-01	Magenta: Medium Limited Range Retention Insert - 4 Pack
30029-01	Gold: High Limited Range Retention Insert - 4 Pack
30012-01	Black: Processing Insert - 4 Pack
08514	Block-Out Spacer - 20 Pack
30013-01	Denture Attachment Housing - 4 Pack
30014-01	Abutment Analog (3.35mm) - 4 Pack



LOCATOR R-Tx® PROCESSING COMPONENTS (CONTINUED)	
PART#	DESCRIPTION
30015-01	Abutment Analog (4.0mm) - 4 Pack
30016-01	Abutment Analog (5.0mm) - 4 Pack
30017-01	Impression Coping - 4 Pack
30018-01	Processing Spacer - 4 Pack
30030-01	Mixed Retention Limited Range Processing Pack - 2 Pack
30031-01	Limited Range Retention Inserts, Processing Pack - 1 Pack
LOCATOR R-Tx® TOOLS	
PART#	DESCRIPTION
30021-01	Retention Insert Tool (Dual Sided)
08007	.050"/1.25mm Hex Driver with Thumb Knob (Shaft 6.4mm)
08008	.050"/1.25mm Hex, Latchlock Driver (Shaft 12.7mm)
04936	.050"/1.25mm Hex, Torque Wrench Insert Square Driver Connection (Shaft 6.4mm)
LOCATOR R-Tx® PATIENT EDUCATION MODELS	
PART#	DESCRIPTION
09665	LOCATOR R-Tx Patient Education Model (Straight)
09677	LOCATOR R-Tx Patient Education Model (Angled)
30053-01	LOCATOR R-Tx 4x Abutment Demonstration Model
LOCATOR R-Tx® SCAN BODIES	
PART#	DESCRIPTION
30060-01-2	LOCATOR Scan Body (2 Pack)
30060-01-4	LOCATOR Scan Body (4 Pack)
30060-01-10	LOCATOR Scan Body (10 Pack)

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-Tx
LOCATOR CHAIRSIDE
MATERIALS & EQUIPMENT

## LOCATOR R-Tx® ABUTMENTS ARE IMPLANT SPECIFIC



TO ORDER, YOU WILL NEED THE FOLLOWING INFORMATION:

1. Name of the implant system
2. Implant diameter
3. Connection type and platform diameter
4. Tissue measurement (at the highest point)

IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

### ACE SURGICAL

#### BONE LEVEL CONNECTION

(Compatible with infinity Octagon Implants)

2.7mm Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
3.3mm Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06

#### TISSUE LEVEL CONNECTION

(Compatible with infinity Octagon Implants)

4.8mm Regular Platform (RP)		31302-01	31302-02	31302-03	31302-04	31302-05	31302-06
-----------------------------	--	----------	----------	----------	----------	----------	----------

#### INTERNAL HEX CONNECTION

(Compatible with infinity Internal HEX Implants)

3.5mm	30200-00 (0.5)*	30200-01*	30200-02*	30200-03*	30200-04*	30200-05*	30200-06*
4.5mm	30201-00 (0.5)*	30201-01*	30201-02*	30201-03*	30201-04*	30201-05*	30201-06*

#### TRI-CHANNEL CONNECTION

(Compatible with infinity TRI-CAM Implants)

3.5mm	30500-00 (0.25)*	30500-01*	30500-02*	30500-03*	30500-04*	30500-05*	30500-06*
4.3mm	30501-00 (0.5)*	30501-01*	30501-02*	30501-03*	30501-04*	30501-05*	30501-06*
5.0mm		30502-01*	30502-02*	30502-03*	30502-04*	30502-05*	30502-06*

#### EXTERNAL HEX CONNECTION

(Compatible with infinity EXTERNAL HEX Implants)

3.5mm		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm (Platform Switch Only)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06

### AVINENT

#### INTERNAL HEX CONNECTION

(Compatible with Biomimetic Ocean and Coral)

3.5mm I.C (Ocean)	32500-00 (0.5)	32500-01	32500-02	32500-03	32500-04	32500-05	32500-06
4.1mm I.C (Ocean & Coral)	32500-00 (0.5)	32500-01	32500-02	32500-03	32500-04	32500-05	32500-06

#### EXTERNAL HEX CONNECTION

(Compatible with Biomimetic Ocean and Coral)

3.5mm E.C		32505-01	32505-02	32505-03	32505-04	32505-05	32505-06
-----------	--	----------	----------	----------	----------	----------	----------

IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## BIOHORIZONS

### INTERNAL HEX CONNECTION

(Compatible with Tapered Plus, Tapered Internal and Single-Stage)

3.0mm		30600-01*	30600-02*	30600-03*	30600-04*	30600-05*	30600-06*
3.5mm	30200-00 (0.5)*	30200-01*	30200-02*	30200-03*	30200-04*	30200-05*	30200-06*
4.5mm	30201-00 (0.5)*	30201-01*	30201-02*	30201-03*	30201-04*	30201-05*	30201-06*
5.7mm		30202-01*	30202-02*	30202-03*	30202-04*	30202-05*	30202-06*

### TAPERED PRO CONICAL CONNECTION

(Compatible with CONELOG)



3.3 mm Narrow Platform		30805-01	30805-02	30805-03	30805-04	30805-05	30805-06
3.8 mm Regular Platform		30806-01	30806-02	30806-03	30806-04	30806-05	30806-06

## BIOMET 3i NOTE: BIOMET 3i Abutments feature a .048 hex drive

### INTERNAL HEX CERTAIN® CONNECTION

(Compatible with 3i T3®, NanoTite™ and OSseotite® Implants)

3.4mm		30300-01*	30300-02*	30300-03*	30300-04*	30300-05*	30300-06*
4.1mm		30301-01*	30301-02*	30301-03*	30301-04*	30301-05*	30301-06*

### EXTERNAL HEX CONNECTION

(Compatible with 3i T3®, NanoTite™ and OSseotite® Implants)

3.4mm		30305-01	30305-02	30305-03	30305-04	30305-05	30305-06
4.1mm		30306-01 (1.2)	30306-02	30306-03	30306-04	30306-05	30306-06

## BLUE SKY BIO

### BONE LEVEL CONNECTION

(Compatible with Quattro Implants)

2.8mm Narrow Platform (NP)		30900-01	30900-02	30900-03	30900-04	30900-05	30900-06
3.3mm Regular Platform (RP)		30901-01	30901-02	30901-03	30901-04	30901-05	30901-06

### TISSUE LEVEL CONNECTION

(Compatible with One Stage Implants)

4.8mm Regular Platform (RP)		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06
6.5mm Wide Platform (WP)		30907-01	30907-02	30907-03	30907-04	30907-05	30907-06

## DENTSPLY

### CONICAL CONNECTION

(Compatible with Ankylos® C/X Implants)

2.8mm Ankylos Platform		30715-01*	30715-02*	30715-03*	30715-04*	30715-05*	30715-06*
------------------------	--	-----------	-----------	-----------	-----------	-----------	-----------

\*Available in Abutment Only configuration. Add "-AO" to order part without processing components.

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-Tx
LOCATOR CHAIRSIDE
MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

### DENTSPLY (continued)

#### (Compatible with AstraTech Osseospeed EV Implants)

3.6mm		30700-01*	30700-02*	30700-03*	30700-04*	30700-05*	30700-06*
4.2mm		30701-01*	30701-02*	30701-03*	30701-04*	30701-05*	30701-06*
4.8mm		30702-01*	30702-02*	30702-03*	30702-04*	30702-05*	30702-06*

#### CONICAL CONNECTION

#### (Compatible with AstraTech Osseospeed TX Implants)

2.9mm Small Platform (Aqua)		30706-01*	30706-02*	30706-03*	30706-04*	30706-05*	30706-06*
3.9mm Large Platform (Lilac)		30707-01	30707-02	30707-03	30707-04	30707-05	30707-06

#### INTERNAL CONICAL CONNECTION

#### (Compatible with OmniTaper EV, Astra Tech EV - OsseoSpeed®, and PrimeTaper EV)

3.6mm		30700-01	30700-02	30700-03	30700-04	30700-05	30700-06
4.2 mm		30701-01	30701-02	30701-03	30701-04	30701-05	30701-06
4.8 mm		30702-01	30702-02	30702-03	30702-04	30702-05	30702-06

#### INTERNAL HEX CONNECTION

#### (Compatible with XiVE® Frialit-2 Implants)

3.4mm		30711-01	30711-02	30711-03	30711-04	30711-05	30711-06
3.8mm		30712-01	30712-02	30712-03	30712-04	30712-05	30712-06
4.5mm		30713-01	30713-02	30713-03	30713-04	30713-05	30713-06
5.5mm		30714-01	30714-02	30714-03	30714-04	30714-05	30714-06

### DIO

#### INTERNAL CONNECTION

#### (Compatible with DIO UF (II))

3.0/3.3mm Narrow		31405-01	31405-02	31405-03	31405-04	31405-05	31405-06
3.8/4.0/4.5/5.0/5.5mm Regular/Wide		31406-01	31406-02	31406-03	31406-04	31406-05	31406-06

### FAIR IMPLANT

#### CONICAL CONNECTION

#### (Compatible with FairTwo)

2.8mm		31900-01	31900-02	31900-03	31900-04	31900-05	31900-06
3.5mm		31901-01	31901-02	31901-03	31901-04	31901-05	31901-06

### HIOSSEN/OSSTEM

#### EXTERNAL HEX CONNECTION

#### (Compatible with US)

3.5mm Mini Platform (MP) (Yellow)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Regular Platform (RP) (Green)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.1mm Wide Platform (WP) (Blue)		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06

IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	DESCRIPTION	CUFF HEIGHT					
		0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm

### HIOSSEN/OSSTEM (continued)

#### CONICAL/MORSE TAPE CONNECTION

(Compatible with ET and TS)

3.35mm Mini Regular Platform (RP) (Yellow)		31601-01	31601-02	31601-03	31601-04	31601-05	31601-06
--------------------------------------------	--	----------	----------	----------	----------	----------	----------

#### INTERNAL HEX/MORSE TAPER CONNECTION

(Compatible with ET and TS)

2.8mm Mini Platform (MP) (Yellow)		31600-01	31600-02	31600-03	31600-04	31600-05	31600-06
-----------------------------------	--	----------	----------	----------	----------	----------	----------

#### INTERNAL OCTAGON/MORSE TAPER CONNECTION

(Compatible with SS)

4.8mm Regular Platform (RP) (Green)		31606-01	31606-02	31606-03	31606-04	31606-05	31606-06
-------------------------------------	--	----------	----------	----------	----------	----------	----------

### IDS

#### INTERNAL HEX CONNECTION (Compatible with MegaGen AnyRidge)

3.5mm Small Platform (SP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Regular Platform (RP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Wide Platform (WP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
5.0mm Super Wide Platform		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06

#### INTERNAL HEX CONNECTION

(Compatible with Legacy1, Legacy2, Legacy3 and Legacy4 Implants)

3.0mm (Blue) (Legacy 2,3,4 Only)		31500-01 (0.75)*	31500-02*	31500-03*	31500-04*	31500-05*	31500-06*
3.5mm (Green)	30200-00 (0.5)*	30200-01*	30200-02*	30200-03*	30200-04*	30200-05*	30200-06*
4.5mm (Purple)	30201-00 (0.5)*	30201-01*	30201-02*	30201-03*	30201-04*	30201-05*	30201-06*
5.7mm (Yellow)		30202-01*	30202-02*	30202-03*	30202-04*	30202-05*	30202-06*

#### TRI-LOBE CONNECTION

(Compatible with RePlant and ReActive Implants)

3.5mm (Purple)	30500-00 (0.25)*	30500-01*	30500-02*	30500-03*	30500-04*	30500-05*	30500-06*
4.3mm (Yellow)	30501-00 (0.5)*	30501-01*	30501-02*	30501-03*	30501-04*	30501-05*	30501-06*
5.0mm (Blue)		30502-01*	30502-02*	30502-03*	30502-04*	30502-05*	30502-06*

#### CONICAL CONNECTION

(Compatible with Interactive Implants)

3.0mm (Purple)		30506-01*	30506-02*	30506-03*	30506-04*	30506-05*	30506-06*
3.4mm (Yellow)		30507-01*	30507-02*	30507-03*	30507-04*	30507-05*	30507-06*

#### INTERNAL OCTAGON CONNECTION

(Compatible with Swish Plus and Swish Tapered Implants)

3.7mm		31500-01*	31500-02*	31500-03*	31500-04*	31500-05*	31500-06*
4.8mm		30906-01	30906-02	30906-03	30906-04	30906-05	30906-06

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-Tx

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## KEYSTONE DENTAL

### CONICAL CONNECTION

(Compatible with Genesis Active)

2.9 mm KGD		31203-01-SB	31203-02-SB	31203-03-SB	31203-04-SB	31203-05-SB	31203-06-SB
------------	--	-------------	-------------	-------------	-------------	-------------	-------------

### TiLOBE CONNECTION

(Compatible with Genesis and PrimaConnex)

3.5mm (Green)		31200-01	31200-02	31200-03	31200-04	31200-05	31200-06
---------------	--	----------	----------	----------	----------	----------	----------

### TiLOBE CONNECTION

(Compatible with PrimaConnex)

4.1mm (Blue)		31201-01	31201-02	31201-03	31201-04	31201-05	31201-06
--------------	--	----------	----------	----------	----------	----------	----------

(Compatible with PrimaConnex and TiLOBEMAXX)

5.0mm (Purple)		31202-01	31202-02	31202-03	31202-04	31202-05	31202-06
----------------	--	----------	----------	----------	----------	----------	----------

## MEGAGEN/IDS

### INTERNAL HEX CONNECTION

(Compatible with MegaGen AnyRidge)

3.5mm Small Platform (SP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Regular Platform (RP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
3.5mm Wide Platform (WP)		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06
5.0mm Super Wide Platform		32201-01	32201-02	32201-03	32201-04	32201-05	32201-06

## MIS

### CONICAL CONNECTION

(Compatible with C1 and V3 Implants)

2.75mm Narrow Platform (NP)		32601-01	32601-02	32601-03	32601-04	32601-05	32601-06
3.15mm Standard Platform (SP)		32602-01	32602-02	32602-03	32602-04	32602-05	32602-06

### INTERNAL HEX CONNECTION (Compatible with SEVEN and M4 Implants)

3.3mm Narrow Platform (NP)		32605-01	32605-02	32605-03	32605-04	32605-05	32605-06
3.75mm Standard Platform (SP)	32606-00 (0.25)	32606-01	32606-02	32606-03	32606-04	32606-05	32606-06
4.5mm Wide Platform (WP)	32607-00 (0.5)	32607-01	32607-02	32607-03	32607-04	32607-05	32607-06

## MOZO GRAU

### EXTERNAL HEX CONNECTION (Compatible with Osseous)

4.1mm		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.0mm		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06

### INTERNAL HEX CONNECTION (Compatible with Inhex)

2.8mm		30706-01*	30706-02*	30706-03*	30706-04*	30706-05*	30706-06*
-------	--	-----------	-----------	-----------	-----------	-----------	-----------

\*Available in Abutment Only configuration. Add "-AO" to order part without processing components.

IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## NEODENT

INTERNAL HEX CONNECTION (Compatible with CM Implants)							
2.5mm		32301-01*	32301-02*	32301-03*	32301-04*	32301-05*	32301-06*
CONICAL/GRAND MORSE CONNECTION (Compatible with GM Implants)							
3.0mm		32700-01*	32700-02*	32700-03*	32700-04*	32700-05*	32700-06*

## NOBEL BIOCARE

CONICAL CONNECTION (Compatible with NobelActive®, NobelParallel™ CC and NobelReplace® CC Implants)							
3.5mm Narrow Platform (NP)		30506-01*	30506-02*	30506-03*	30506-04*	30506-05*	30506-06*
4.3mm Regular Platform (RP)		30507-01*	30507-02*	30507-03*	30507-04*	30507-05*	30507-06*
5.0mm Wide Platform (WP)		30508-01*	30508-02*	30508-03*	30508-04*	30508-05*	30508-06*

EXTERNAL HEX CONNECTION (Compatible with Brånemark and NobelSpeedy® Groovy Implants)							
3.5mm Narrow Platform (NP)		30510-01	30510-02	30510-03	30510-04	30510-05	30510-06
4.1mm Regular Platform (RP)		30511-01	30511-02	30511-03	30511-04	30511-05	30511-06
5.1mm Wide Platform (WP)		30512-01	30512-02	30512-03	30512-04	30512-05	30512-06

TRI-CHANNEL CONNECTION (Compatible with NobelReplace®, ReplaceSelect™ and NobelSpeedy® Replace Implants)							
3.5mm Narrow Platform (NP)	30500-00 (0.25)*	30500-01*	30500-02*	30500-03*	30500-04*	30500-05*	30500-06*
4.0mm/4.3mm Regular Platform (RP)	30501-00 (0.5)*	30501-01*	30501-02*	30501-03*	30501-04*	30501-05*	30501-06*
5.0mm Wide Platform (WP)		30502-01*	30502-02*	30502-03*	30502-04*	30502-05*	30502-06*

## OCO BIOMEDICAL

INTERNAL HEX CONNECTION (Compatible with Engage® Implants)							
3.0mm		31500-01 (0.75)*	31500-02*	31500-03*	31500-04*	31500-05*	31500-06*
3.5mm	30200-00 (0.5)*	30200-01*	30200-02*	30200-03*	30200-04*	30200-05*	30200-06*
4.5mm	30201-00 (0.5)*	30201-01*	30201-02*	30201-03*	30201-04*	30201-05*	30201-06*

EXTERNAL HEX CONNECTION (Compatible with Engage® Implants)							
5.7mm		30202-01*	30202-02*	30202-03*	30202-04*	30202-05*	30202-06*

## OSTEOREADY/LITTLE IMPLANT COMPANY

INTERNAL HEX CONNECTION (Compatible with Marc Nevins®, Pamela Ray™ and Performance®)							
3.5mm Diameter Platform	32100-00 (0.5)	32100-01	32100-02	32100-03	32100-04	32100-05	32100-06

LOCATOR Root  
LOCATOR Implants  
LOCATOR Abutments  
LOCATOR Bar  
LOCATOR FIXED  
REGAL HYGIENE  
LOCATOR R-Tx  
LOCATOR CHAIRSIDE  
MATERIALS & EQUIPMENT

IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
DESCRIPTION	CUFF HEIGHT						
	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## SIC

### INTERNAL HEX CONNECTION

(Compatible with SICace/SICmax)

3.3mm		31805-01	31805-02	31805-03	31805-04	31805-05	31805-06
4.2 mm		31806-01	31806-02	31806-03	31806-04	31806-05	31806-06

### CONICAL CONNECTION

(Compatible with SICvantage max Implants)

2.2mm		31800-01 (1.5mm)	31800-02	31800-03	31800-04	31800-05	31800-06
-------	--	---------------------	----------	----------	----------	----------	----------

## SOUTHERN IMPLANTS

### TRI-NEX CONNECTION

(Compatible with Tri-Nex™ Implants)

3.5mm Tri-Nex Platform	30500-00 (0.25mm)	30500-01	30500-02	30500-03	30500-04	30500-05	30500-06
4.3mm Tri-Nex Platform	30501-00 (0.5mm)	30501-01	30501-02	30501-03	30501-04	30501-05	30501-06
5.0mm Tri-Nex Platform		30502-01	30502-02	30502-03	30502-04	30502-05	30502-06

## STRAUMANN

### BONE LEVEL

(Compatible with Roxolid® SLA® and Roxolid® SLActive® Implants)

3.3mm Narrow Crossfit (NC)		30401-01*	30401-02*	30401-03*	30401-04*	30401-05*	30401-06*
4.1mm/4.8mm Regular Crossfit (RC)		30402-01*	30402-02*	30402-03*	30402-04*	30402-05*	30402-06*

(Compatible with BLX Implants)

Regular Base (RB)		30410-01*	30410-02*	30410-03*	30410-04*	30410-05*	30410-06*
Wide Base (WB)		30410-01*	30410-02*	30410-03*	30410-04*	30410-05*	30410-06*

### TISSUE LEVEL

(Compatible with Roxolid SLA and Roxolid SLActive Implants)

3.5mm Narrow Neck Crossfit (NNC)		30405-01*	30405-02*	30405-03*	30405-04*	30405-05*	30405-06*
4.1mm Regular Neck (RN)		30406-01*	30406-02*	30406-03*	30406-04*	30406-05*	30406-06*
4.8mm Wide Neck (WN)		30407-01* (1.5)	30407-02*	30407-03*	30407-04*	30407-05*	30407-06*

## SOUTHERN IMPLANTS

### EXTERNAL HEX CONNECTION

4.05 Diameter Platform		31703-01	31703-02	31703-03	31703-04	31703-05	31703-06
5.0 Diameter Platform		31704-01	31704-02	31704-03	31704-04	31704-05	31704-06
6.0 Diameter Platform		31705-01	31705-02	31705-03	31705-04	31705-05	31705-06

\*Available in Abutment Only configuration. Add "-AO" to order part without processing components.



IMPLANT COMPANY	LOCATOR R-Tx® ATTACHMENT SYSTEM (Includes Abutment and Processing Components)						
	CUFF HEIGHT						
DESCRIPTION	0.0mm	1.0mm	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm

## SWEDEN & MARTINA

### INTERNAL HEX CONNECTION

(Compatible with Premium/Prama)

3.3mm		31110-01	31110-02	31110-03	31110-04	31110-05	31110-06
3.8mm		31111-01	31111-02	31111-03	31111-04	31111-05	31111-06
4.25mm		31112-01	31112-02	31112-03	31112-04	31112-05	31112-06
5.0mm		31113-01	31113-02	31113-03	31113-04	31113-05	31113-06

(Compatible with Syra)

4.1mm		31114-01	31114-02	31114-03	31114-04	31114-05	31114-06
-------	--	----------	----------	----------	----------	----------	----------

### EXTERNAL HEX CONNECTION

(Compatible with Outlink)

3.3mm		31105-01	31105-02	31105-03	31105-04	31105-05	31105-06
4.1mm		31106-01	31106-02	31106-03	31106-04	31106-05	31106-06

## ZIMMER

### INTERNAL HEX CONNECTION

(Compatible with Tapered Screw-Vent®, Screw-Vent®, Trabecular Metal™ Implants and AdVent® 5.7mm Platform Implant)

3.5mm	30200-00 (0.5)*	30200-01*	30200-02*	30200-03*	30200-04*	30200-05*	30200-06*
4.5mm	30201-00 (0.5)*	30201-01*	30201-02*	30201-03*	30201-04*	30201-05*	30201-06*
5.7mm		30202-01*	30202-02*	30202-03*	30202-04*	30202-05*	30202-06*

### SPLINE CONNECTION

(Compatible with Spline Reliance and Cylinder Implants)

3.25mm		30205-01	30205-02	30205-03	30205-04	30205-05	30205-06
3.75mm		30206-01	30206-02	30206-03	30206-04	30206-05	30206-06
4.0mm		30206-01	30206-02	30206-03	30206-04	30206-05	30206-06
5.0mm		30207-01	30207-02	30207-03	30207-04	30207-05	30207-06

LOCATOR Root  
LOCATOR Implants  
LOCATOR Abutments  
LOCATOR Bar  
LOCATOR FIXED  
REGAL HYGIENE  
LOCATOR R-Tx  
LOCATOR CHAIRSIDE  
MATERIALS & EQUIPMENT

## Full-Arch Materials and Tools

# LOCATOR CHAIRSIDE®

The perfect set of materials and tools to support LOCATOR and full-arch procedures chairside.



### LOCATOR CHAIRSIDE® IMPRESSION MATERIALS

PART#	DESCRIPTION
09614	LOCATOR CHAIRSIDE Impression Material, Light Body, 2-50ml cartridges, 6 mixing tips (MAGENTA)
09615	LOCATOR CHAIRSIDE Impression Material, Medium Body, 2-50ml cartridges, 6 mixing tips (DARK BLUE)
09616	LOCATOR CHAIRSIDE Impression Material, Heavy Body, 2-50ml cartridges, 6 mixing tips (YELLOW)

### LOCATOR CHAIRSIDE® BITE REGISTRATION MATERIAL

PART#	DESCRIPTION
09572	LOCATOR CHAIRSIDE Bite Registration Material, 2-50ml cartridges, 6 mixing tips

### LOCATOR CHAIRSIDE® ATTACHMENT PROCESSING MATERIAL

PART#	DESCRIPTION
09566	LOCATOR CHAIRSIDE Attachment Processing Material Normal Set 8g Syringe (Includes: 15 Brown Mixing Tips and 15 Angled Tips)

**LOCATOR CHAIRSIDE® ATTACHMENT PROCESSING MATERIAL (continued)**

PART#	DESCRIPTION
90-00130	LOCATOR CHAIRSIDE Attachment Processing Material Fast Set 8g Syringe (Includes: 15 Brown Mixing Tips and 15 Angled Tips)
09565	LOCATOR CHAIRSIDE Attachment Processing Material Normal Set 18mL Cartridge (Includes: 10 Blue Mixing Tips)
90-00132	LOCATOR CHAIRSIDE Attachment Processing Material Fast Set 18mL Cartridge (Includes: 10 Blue Mixing Tips)
09125-15	LOCATOR CHAIRSIDE Attachment Processing Material Brown Mixing Tips (15) (To use with 8g syringe)
09973-15	LOCATOR CHAIRSIDE Attachment Processing Material Blue Mixing Tips (15) (To use with 18mL cartridge)
91428-60-M	LOCATOR CHAIRSIDE APM intraoral tip

**LOCATOR CHAIRSIDE® TOOLS AND PREP AND POLISH KIT**

PART#	DESCRIPTION
09582	LOCATOR CHAIRSIDE Denture Prep and Polish Kit (Includes: Recess, Trim, Undercut, Grind and Vent Burs and Polisher)
09576	LOCATOR CHAIRSIDE Recess Bur
09577	LOCATOR CHAIRSIDE Undercut Bur
09578	LOCATOR CHAIRSIDE Vent Bur
09579	LOCATOR CHAIRSIDE Trim Bur
09583	LOCATOR CHAIRSIDE Grind Bur
09580	LOCATOR CHAIRSIDE Polishing Bur (Includes: 5 Polishers)
09702	LOCATOR CHAIRSIDE Trepine Bur (Not included with the Denture Prep and Polish Kit)

**LOCATOR CHAIRSIDE® SOFT RELINE MATERIAL KIT**

PART#	DESCRIPTION
09573	LOCATOR CHAIRSIDE Soft Reline Material Kit (Includes: 52g Cartridge, 10mL Primer Bottle, 10 Mixing Tips, Plastic Cup, Brush, Spatula, Coarse Polishing Point and Finishing Wheel)

**LOCATOR CHAIRSIDE® BLOCK OUT MATERIAL**

PART#	DESCRIPTION
09574	LOCATOR CHAIRSIDE Block Out Material (Includes: (4) 2mL Syringes and 20 Mixing Tips)

**LOCATOR CHAIRSIDE® DENTURE REMOVAL TOOL**

PART#	DESCRIPTION
09549	LOCATOR CHAIRSIDE Patient Removal Tool (Includes 10 Tools)

LOCATOR Root  
 LOCATOR Implants  
 LOCATOR Abutments  
 LOCATOR Bar  
 LOCATOR FIXED  
 REGAL HYGIENE  
 LOCATOR R-TX  
 LOCATOR CHAIRSIDE  
 MATERIALS & EQUIPMENT



# DENTAL MATERIALS & EQUIPMENT

# TABLE OF CONTENTS

## DENTAL MATERIALS & EQUIPMENT

CARIES DETECTION & DESENSITIZING.....	62
AIR INHIBITORS, BONDING AGENTS AND ETCHANTS.....	62
PRELUDE BONDING AGENT.....	62
ETCHANT.....	63
COMPOSITES.....	64
ZNANO UNIVERSAL FLOW.....	64
ZNANO - UNIVERSAL NANOCOMPOSITE.....	65
BULK EZ PLUS - BULKFILL COMPOSITE.....	66
BULK EZ - BULKFILL COMPOSITE.....	66
ROCK CORE - CORE BUILD-UP.....	67
STARFILL 2B - DUAL-CURE COMPOSITE.....	67
ARIA MICROFILL - FLOWABLE COMPOSITE.....	67
STARFLOW - FLOWABLE COMPOSITE.....	68
STARFLOW PV - FLOWABLE COMPOSITE.....	69
ACCOLADE - HIGH VISCOSITY FLOWABLE COMPOSITE.....	69
ACCOLADE SUPER RADIOPAQUE (SRO) - FLOWABLE COMPOSITE.....	69
ACCOLADE PV VENEER PLACEMENT SYSTEM.....	70
ACCOLADE PV VENEER - FLOWABLE COMPOSITE.....	70
ACCOLADE PV - TRY-IN COMPOSITE PASTE.....	70
TIPS.....	71
ACCOLADE OP MASK.....	71
CEMENT.....	72
CEMEZ - ADHESIVE RESIN CEMENT.....	72
IMPRESSION MATERIAL.....	72
START VPS - REGULAR SET.....	72
FIRST HALF - MEDIUM SET.....	73
FIRST QUARTER - FAST SET.....	74
VPS ALGINATE SUBSTITUTE.....	74
TIPS AND ACCESSORIES.....	74
MATRIX SYSTEMS.....	75
CONTACT MATRIX.....	75
WEDGES.....	76

STAIN REMOVAL.....	77
FLASHBUSTER COMPOSITE BURS.....	77
POSTS.....	77
ICELIGHT - LIGHT TRANSMITTING COMPOSITE POST.....	77
ICEPOST - COMPOSITE POST.....	77
REAMERS.....	78
PROVISIONAL COMPOSITES.....	78
TURBOTEMP EZ.....	78
TURBOTEMP 2.....	79
TURBOTEMP 3.....	79
TIPS AND ACCESSORIES.....	79
TEMPORARY CEMENT.....	79
IVERI WHITENING.....	80
TAKE-HOME KITS.....	80
TAKE-HOME REFILL SYRINGES.....	81
AIR ABRASION.....	82
PREPSTART UNITS AND ACCESSORIES.....	82
PREPSTART NOZZLES.....	82
RETRACTORS.....	83
SANDBLASTING.....	83
MICROETCHER UNITS AND ACCESSORIES.....	83
MICROETCHER NOZZLES.....	84
INSTALLATION KITS.....	84
SILJET.....	85
ABRASIVES.....	85
DUST CONFINEMENT.....	86
DUST CABINETS.....	86
MICROCAB REPLACEMENT PARTS.....	86
MACROCAB REPLACEMENT PARTS.....	86
DC-2 REPLACEMENT PARTS.....	87

## CARIES DETECTION & DESENSITIZING



PART#	DESCRIPTION
80010	Caries Finder Red (10 ml)
80004	Caries Finder Green (10 ml)
87001	MicroPrime B Type (10 ml)
90814	MicroPrime G Type (10 ml)
94263	MicroPrime G Type Unit Dose (0.2 ml), 50 pack
90-00249	Microprime G Gel, (1 X 4 ml syringe)
90-00250	Microprime G Gel, (4 X 1 ml syringe)

## AIR INHIBITORS, BONDING AGENTS AND ETCHANTS



PRELUDE BONDING AGENT	
PART#	DESCRIPTION
91653	Prelude One Universal Adhesive (5 ml)
91657	Prelude One Unit Dose (0.1 ml), 50 pack
90990	Prelude Self-Etch Kit: Primer (5 ml), Adhesive (5 ml), Link (5 ml), 10 Mixing Wells, 40 Applicator Brushes

**PRELUDE BONDING AGENT (CONTINUED)**

PART#	DESCRIPTION
90994	Prelude Total-Etch Kit: Adhesive (5 ml), MicroPrime G (5 ml), SureEtch (5 g), 10 Mixing Wells, 40 Applicator Brushes
91024	Prelude Primer/Adhesive 2 pack (5 ml ea.)
90972	Prelude Primer (5 ml)
90974	Prelude Adhesive (5 ml)
90976	Prelude Dual/Self-Cure Link (5 ml)
94070	Prelude Primer Unit Dose (0.2 ml), 30 pack
94071	Prelude Adhesive Unit Dose (0.2 ml), 30 pack
94072	Prelude Link Unit Dose (0.2 ml), 30 pack
93139	H2Orthobond (5ml)
93410	Z-Bond (5ml)

**AIR INHIBITORS, BONDING AGENTS AND ETCHANTS (CONTINUED)**

PART#	DESCRIPTION
87200	Liquid Lens (5 ml), 12 Tips
90752	Dam Cool (5 ml), 12 Tips
90461	No Stick (10 ml)
89104	S-Bond (10 ml)
89107	E-Bond (10 ml)

**ETCHANT**

PART#	DESCRIPTION
89463	SureEtch Gel (5ml), 12 Tips (20 gauge),
85072	SureEtch Liquid (10 ml)

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

**MATERIALS & EQUIPMENT**

## COMPOSITES

### ZNANO - UNIVERSAL FLOW



PART#	DESCRIPTION
90-00137	ZNano Universal Flowable, High Flow, A1/B1
90-00138	ZNano Universal Flowable, High Flow, A2/B2
90-00139	ZNano Universal Flowable, High Flow, A3/A3.5/B3
90-00140	ZNano Universal Flowable, High Flow, C2/C3
90-00141	ZNano Universal Flowable, High Flow, BW/OP
90-00142	ZNano Universal Flowable, High Flow, A1/B1
90-00143	ZNano Universal Flowable, High Flow, A2/B2
90-00144	ZNano Universal Flowable, High Flow, A3/A3.5/B3
90-00145	ZNano Universal Flowable, High Flow, C2/C3
90-00146	ZNano Universal Flowable, High Flow, BW/OP



## ZNANO - UNIVERSAL NANOCOMPOSITE



PART#	DESCRIPTION
90-00112	ZNano Shade-Adaptive Universal Nanocomposite A1/B1 Syringe Refill (4 gm)
90-00113	ZNano Shade-Adaptive Universal Nanocomposite A2/B2 Syringe Refill (4 gm)
90-00114	ZNano Shade-Adaptive Universal Nanocomposite A3/A3.5/B3 Syringe Refill (4 gm)
90-00115	ZNano Shade-Adaptive Universal Nanocomposite Bleach/Opaque Syringe Refill (4 gm)
90-00116	ZNano Shade-Adaptive Universal Nanocomposite C2/C3 Syringe Refill (4 gm)
90-00118	ZNano Shade-Adaptive Universal Nanocomposite A1/B1 Unit Dose Refill (20x 0.2 gm)
90-00119	ZNano Shade-Adaptive Universal Nanocomposite A2/B2 Unit Dose Refill (20x 0.2 gm)
90-00120	ZNano Shade-Adaptive Universal Nanocomposite A3/A3.5/B3 Unit Dose Refill (20x 0.2 gm)
90-00121	ZNano Shade-Adaptive Universal Nanocomposite Bleach/Opaque Unit Dose Refill (20x 0.2 gm)
90-00122	ZNano Shade-Adaptive Universal Nanocomposite C2/C3 Unit Dose Refill (20x 0.2 gm)
95808	Restorative Materials Trial Kit (Includes samples of Microprime G, Prelude One, and ZNano)



LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## COMPOSITES (CONTINUED)

### BULK EZ PLUS - BULKFILL COMPOSITE



PART#	DESCRIPTION
90-00171	Bulk EZ PLUS A1/B1 Refill, 1x 6 gm syringe, 6 tips (19 gauge)
90-00172	Bulk EZ PLUS A2/B2 Refill, 1x 6 gm syringe, 6 tips (19 gauge)
90-00173	Bulk EZ PLUS A3/A3.5/B3 Refill, 1x 6 gm syringe, 6 tips (19 gauge)
90-00174	Bulk EZ PLUS C2/C3 Refill, 1x 6 gm syringe, 6 tips (19 gauge)
90-00175	Bulk EZ PLUS C2/C3 Refill, 1x 6 gm syringe, 6 tips (19 gauge)
90-00176	Bulk EZ PLUS Core White Refill, 1x 6 gm syringe, 6 tips (19 gauge)
90-00177	Bulk EZ PLUS Trial Kit, 1x 6gm syringe (A1/B1), 1x 6gm syringe (A2/B2), 1x 6gm syringe (A3/A3.5/B3), 1x 6gm syringe (C2C3), 1x 6gm syringe (BL/OP), 1x 6gm syringe (CORE White), 18 Tips (17 gauge), 18 Tips (19 gauge)

### BULK EZ - BULK FILL COMPOSITE



PART#	DESCRIPTION
94778	Bulk EZ Trial Kit 3 Syringes (6 gm ea.) -- A1, A2, A3 -- 20 Automix Tips (10 each: 17 and 19 gauge)
94775	Bulk EZ Refill A1, 2 Syringes (6 gm ea.), 6 Tips (17 gauge), 6 Tips (19 gauge)
94776	Bulk EZ Refill A2, 2 Syringes (6 gm ea.), 6 Tips (17 gauge), 6 Tips (19 gauge)
94777	Bulk EZ Refill A3, 2 Syringes (6 gm ea.), 6 Tips (17 gauge), 6 Tips (19 gauge)
94792	Bulk EZ Tips (17 Gauge), 25 pack
94793	Bulk EZ Tips (19 Gauge), 25 pack

## ROCK CORE - CORE BUILD-UP

PART#	DESCRIPTION
90838	Rock Core, White (25 ml), 6 Mixing Tips w/Intraoral Tip
93831	Rock Core, A-2 (25 ml), 6 Mixing Tips w/Intraoral Tip
90605	Rock Core, Automix Snap-On Intraoral Tip, 100 pack
90130	Rock Core 25 ml Dispenser Gun

## STARTFILL 2B - DUAL-CURE COMPOSITE



PART#	DESCRIPTION
90150	StartFill 2B (10 gm), 6 Tips
91686	StartFill 2B Mixing tip w/Intraoral Tip, 12 pack
91430	StartFill 2B Mixing tip w/Intraoral Tip, 48 pack

## ARIA MICROFILL - FLOWABLE COMPOSITE



PART#	DESCRIPTION
89360	Aria A1 (5 gm), 6 Tips (18 gauge)
89361	Aria A2 (5 gm), 6 Tips (18 gauge)
89362	Aria A3.5 (5 gm), 6 Tips (18 gauge)
89356	Aria A1 (1 gm), 6 Tips (18 gauge)
89357	Aria A2 (1 gm), 6 Tips (18 gauge)
89358	Aria A3.5 (1 gm), 6 Tips (18 gauge)

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## COMPOSITES (CONTINUED)

### STARTFLOW - FLOWABLE COMPOSITE



PART#	DESCRIPTION
85071	Intro Kit, A2, A3.5, Incisal syringes (5 gm), 100 Tips (20 gauge)
89421	A00 OP (5 gm), 6 Tips (20 gauge)
89397	A0 (5 gm), 6 Tips (20 gauge)
85051	A1 (5 gm), 6 Tips (20 gauge)
85052	A2 (5 gm), 6 Tips (20 gauge)
85053	A3 (5 gm), 6 Tips (20 gauge)
85054	A3.5 (5 gm), 6 Tips (20 gauge)
90750	A4 (5 gm), 6 Tips (20 gauge)
85041	A5 (5 gm), 6 Tips (20 gauge)
90471	B1 (5 gm), 6 Tips (20 gauge)
85055	B2 (5 gm), 6 Tips (20 gauge)
85056	C2 (5 gm), 6 Tips (20 gauge)
85057	C4 (5 gm), 6 Tips (20 gauge)
85058	Dentin (5 gm), 6 Tips (20 gauge)
85059	Incisal (5 gm), 6 Tips (20 gauge)
89465	A00 OP (1 gm), 6 Tips (20 gauge)
89467	A1 (1 gm), 6 Tips (20 gauge)
89468	A2 (1 gm), 6 Tips (20 gauge)
89469	A3 (1 gm), 6 Tips (20 gauge)
89470	A3.5 (1 gm), 6 Tips (20 gauge)

## STARTFLOW PV - FLOWABLE COMPOSITE

PART#	DESCRIPTION
89035	Translucent (5 gm), 6 Tips (20 gauge)
89037	Light (3 gm), 6 Tips (20 gauge)
89039	Extra Light (3 gm), 6 Tips (20 gauge)
89041	White Opaquer (3 gm), 6 Tips (20 gauge)

## ACCOLADE - HIGH VISCOSITY FLOWABLE COMPOSITE



PART#	DESCRIPTION
90553	A00 OP (5 gm), 6 Tips (18 gauge)
90546	A1 (5 gm), 6 Tips (18 gauge)
90547	A2 (5 gm), 6 Tips (18 gauge)
90548	A3 (5 gm), 6 Tips (18 gauge)
90549	A3.5 (5 gm), 6 Tips (18 gauge)
90840	A4 (5 gm), 6 Tips (18 gauge)
90550	B1 (5 gm), 6 Tips (18 gauge)

## ACCOLADE SUPER RADIOPAQUE (SRO) - FLOWABLE COMPOSITE

PART#	DESCRIPTION
91415	A2 (1 gm), 6 Tips (18 gauge)
91417	Extra Light (1 gm), 6 Tips (18 gauge)
93829	White Opaquer (1 gm), 6 Tips (18 gauge)
91387	A2 (5 gm), 6 Tips (18 gauge)
91389	Extra Light (5 gm), 6 Tips (18 gauge)
91388	A5 (5 gm), 6 Tips (18 gauge)
93830	White Opaquer (5 gm), 6 Tips (18 gauge)

LOCATOR Root

LOCATOR  
Implants

LOCATOR  
Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL  
HYGIENE

LOCATOR R-TX

LOCATOR  
CHAIRSIDE

MATERIALS &  
EQUIPMENT

**COMPOSITES (CONTINUED)****ACCOLADE PV VENEER PLACEMENT SYSTEM**

<b>PART#</b>	<b>DESCRIPTION</b>
91216	Accolade PV Kit with Try-In Pastes PV Composite & Try-in Composite Syringes- Translucent (5 gm), Light (3 gm), Extra Light (3 gm), White Opaquer (3 gm), Brown (3 gm), Yellow (3 gm), Prelude Adhesive (5 ml), S-Bond (5 ml)

**ACCOLADE PV VENEER - FLOWABLE COMPOSITE**

<b>PART#</b>	<b>DESCRIPTION</b>
90560	Translucent (5 gm), 6 Tips (18 gauge)
90561	Light (3 gm), 6 Tips (18 gauge)
90559	Extra Light (3 gm), 6 Tips (18 gauge)
90562	White Opaquer (3 gm), 6 Tips (18 gauge)
91231	Brown (3 gm), 6 Tips (18 gauge)
91232	Yellow (3 gm), 6 Tips (18 gauge)

**ACCOLADE PV - TRY-IN COMPOSITE PASTE**

<b>PART#</b>	<b>DESCRIPTION</b>
91233	Translucent (5 gm), 6 Tips (18 gauge)
91234	Light (3 gm), 6 Tips (18 gauge)
91235	Extra Light (3 gm), 6 Tips (18 gauge)
91236	White Opaquer (3 gm), 6 Tips (18 gauge)
91237	Brown (3 gm), 6 Tips (18 gauge)
91238	Yellow (3 gm), 6 Tips (18 gauge)

## TIPS



PART#	DESCRIPTION
87211	Pink Tips (18 gauge), 20 pack
87201	Pink Tips (18 gauge), 50 pack
87206	Pink Tips (18 gauge), 100 pack
87212	Orange Tips (20 gauge), 20 pack
87202	Orange Tips (20 gauge), 50 pack
87204	Orange Tips (20 gauge), 100 pack
91428-60	LOCATOR CHAIRSIDE APM intraoral tip
90-00251	Black Flocked Tip, 22ga Refill 20pk

## ACCOLADE OP MASK

PART#	DESCRIPTION
93597	Tooth Color (2 gm), 6 Tips (25 gauge)
93598	White (2 gm), 6 Tips (25 gauge)

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## CEMENT



### CEMEZ - ADHESIVE RESIN CEMENT

PART#	DESCRIPTION
90-00010	Cem EZ Starter Kit: 3 Syringes (5 gm ea.), Translucent, Warm, White Opaque, 24 Automix Tips, Sure Etch Gel (5 ml), 12 Etch Tips, Prelude One (5 ml), S-Bond (5 ml), 40 Applicator Brushes, 1 Reusable Well
90-00011	Refill, Translucent, 2 Syringes (5 gm ea.), 16 Automix Tips
90-00012	Refill, Warm, 2 Syringes (5 gm ea.), 16 Automix Tips
90-00013	Refill, White Opaque, 2 Syringes (5 gm ea.), 16 Automix Tips
90-00025	Cement Automix Tips (24)

## IMPRESSION MATERIAL



### START VPS - REGULAR SET

PART#	DESCRIPTION
80061-01	Ultra Light Econo (50 ml Cartridges), 24 pack
80009-01	Light Body Econo (50 ml Cartridges), 24 pack
80007-01	Heavy Flex Econo (50 ml Cartridges), 24 pack



## START VPS - REGULAR SET (CONTINUED)

PART#	DESCRIPTION
80006-01	Heavy Bite Econo (50 ml Cartridges), 24 pack
80065-01	Ultra Heavy Econo (50 ml Cartridges), 24 pack
80037-01	Clear Bite Econo (50 ml Cartridges), 24 pack
93718	Vanilla Bite Econo (50 ml Cartridges), 24 pack
93354	AlphaBite Econo (50 ml Cartridges), 24 pack
80060-01	Ultra Light (50 ml Cartridges), 4 pack, 12 Automix Tips
80011-01	Light Body (50 ml Cartridges), 4 pack, 12 Automix Tips
80013-01	Heavy Flex (50 ml Cartridges), 4 pack, 12 Automix Tips
80014-01	Heavy Bite (50 ml Cartridges), 4 pack, 12 Automix Tips
80062-01	Ultra Heavy (50 ml Cartridges), 4 pack, 12 Automix Tips
80015-01	Clear Bite (50 ml Cartridges), 4 pack, 12 Automix Tips
80024-01	Heavy Bite (50 ml Cartridges), 2 pack, 6 Automix Tips
80025-01	Clear Bite (50 ml Cartridges), 2 pack, 6 Automix Tips
93717	Vanilla Bite (50 ml Cartridges), 2 pack, 6 Automix Tips
93352	AlphaBite (50 ml Cartridges), 2 pack, 6 Automix Tips
80016	Putty (Base 400 gm & Catalyst 400 gm)

## FIRST HALF - MEDIUM SET

PART#	DESCRIPTION
90650	Heavy Body Lime Green Econo Kit, (50 ml Cartridges), 24 pack
90469	Light Body Black Econo Kit, (50 ml Cartridges), 24 pack
91243	Light Body Orange, (50 ml Cartridges), 24 pack
91124	Monophase Ortho Econo Kit, (50 ml Cartridges), 24 pack
90649	Heavy Body Lime Green (50 ml Cartridges), 4 pack, 12 Automix Tips
90468	Light Body Black (50 ml Cartridges), 4 pack, 12 Automix Tips
91242	Light Body Orange (50 ml Cartridges), 4 pack, 12 Automix Tips
91123	Monophase Ortho (50 ml Cartridges), 4 pack, 12 Automix Tips

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

**FIRST HALF - MEDIUM SET (CONTINUED)**

<b>PART#</b>	<b>DESCRIPTION</b>
93133	Putty (Base 400 gm & Catalyst 400 gm)

**FIRST QUARTER - FAST SET**

<b>PART#</b>	<b>DESCRIPTION</b>
89462	Light Body Econo (50 ml Cartridges), 24 pack
90810	Heavy Body Econo (50 ml Cartridges), 24 pack
80008-01	Monophase Econo (50 ml Cartridges), 24 pack
90902	Chocolate Mint Monophase Econo (50 ml Cartridges), 24 pack
89272	Light Body (50 ml Cartridges), 4 pack, 12 Automix Tips
90809	Heavy Body (50 ml Cartridges), 4 pack, 12 Automix Tips
89387	Monophase (50 ml Cartridges), 4 pack, 12 Automix Tips
90900	Chocolate Mint Monophase (50 ml Cartridges), 4 pack, 12 Automix Tips

**VPS ALGINATE SUBSTITUTE**

<b>PART#</b>	<b>DESCRIPTION</b>
93553	Noginate (50 ml Cartridges), 2 pack, 6 Automix Tips
93554	Noginate Econo (50 ml Cartridges), 24 pack

**TIPS AND ACCESSORIES**

<b>PART#</b>	<b>DESCRIPTION</b>
89331-00	Dispenser Gun
93112	Mixing Tip w/Yellow Intraoral Tip, 12 pack (for Light Body Impression Material and Rock Core)
93113	Mixing Tip w/Yellow Intraoral Tip, 48 pack (for Light Body Impression Material and Rock Core)
90002	Mixing Tips, 12 pack (for Medium and Heavy Body Impression Material)
90003	Mixing Tips, 48 pack (for Medium and Heavy Body Impression Material)
91667	Tray Link Adhesive (10 ml)
93735	Mojo II Syringes, 100 pack

## MATRIX SYSTEMS



### CONTACT MATRIX

PART#	DESCRIPTION
90066	Contact Matrix Trial Kit: 2 Contact Rings (1 Inward, 1 Outward), 100 Thin Flex Matrices (50 each: small and large)
89388	Contact Matrix Clinical Kit: 6 Contact Rings (3 Inward, 3 Outward), 200 Thin Flex Matrices (100 each: small and large)
89439	Contact Matrix Clinical Kit Plus: Contact Matrix Clinical Kit, Contact Ring Pliers, Mega Grip Matrix Forceps
91474	Mega Ring Trial Kit: 2 Mega Rings (1 Inward, 1 Outward), 100 Thin Flex Matrices (50 each: small and large)
91475	Mega Ring Clinical Kit: 6 Mega Rings (3 Inward, 3 Outward), 200 Thin Flex Matrices (100 each: small and large)
91476	Mega Ring Clinical Kit Plus: Mega Ring Clinical Kit, Mega Ring Pliers & Mega Grip Matrix Forceps
88038	Contact Ring Refill - 2 rings (1 Inward, 1 Outward)
89419	Contact Ring Inward Refill, 2 rings
89507	Contact Ring Outward Refill, 2 rings
91303	Mega Ring Inward Refill, 2 rings
91304	Mega Ring Outward Refill, 2 rings
91305	Mega Ring Refill - 2 rings (1 Inward, 1 Outward)
94271	Mega V Ring Clinical Kit Plus: Mega V Clinical Kit, Matrix Forceps & Contact Ring Pliers
94272	Mega V Ring Clinical Kit: 4 Mega V Rings, 75 Mega Wedges (25 each: small, medium, large), 200 Thin Flex Matrices (100 each: small and large)

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

**CONTACT MATRIX (CONTINUED)**

<b>PART#</b>	<b>DESCRIPTION</b>
94270	Mega V Ring Clinical Trial Kit: 2 Mega V Rings, 75 Mega Wedges (25 each: small, medium, large), 100 Thin Flex Matrices (50 each: small and large)
94269	Mega V Ring Refill, 2 Rings
91074	Ultra Thin Flex Matrices, Large, 100 pack
91139	Ultra Thin Flex Matrices, Small, 100 pack
89434	Thin Flex Matrices, Large, 100 pack
89432	Thin Flex Matrices, Small, 100 pack
90071	Stiff Flex Matrices, Large, 100 pack
90070	Stiff Flex Matrices, Small, 100 pack
90679	Subgingival Matrices, 50 pack
90777	Ring Pliers Plus
90486	Mega Grip Forceps

**WEDGES**

<b>PART#</b>	<b>DESCRIPTION</b>
91015	ContactWEDGE Kit: 85 Wedges (small), 85 Wedges (medium), 85 Wedges (large)
91016	ContactWEDGE Refill Small, 85 pack
91017	ContactWEDGE Refill Medium, 85 pack
91018	ContactWEDGE Refill Large, 85 pack
91298	ContactWEDGE Pliers
94810	Mega Wedges, Large, Green, 100 pack
94808	Mega Wedges, Medium, Orange, 100 pack
94806	Mega Wedges, Small, Blue, 100 pack

## STAIN REMOVAL

### FLASHBUSTER COMPOSITE BURS

PART#	DESCRIPTION
2540	Intro Kit 6 Burs: Pointed Cone, Flame, Cyl. Pointed, Extended Flame, Long Flame, Bullet
2501	Pointed Cone - Standard 6 Burs
2502	Flame - Standard 6 Burs
2503	Rev. Cone - Standard 6 Burs
2504	Cyl. Pointed - Standard 6 Burs
2505	Extended Flame - Paro 6 Burs
2506	Long Flame - Paro 6 Burs
4001	Rev. Cone - Jumbo 6 Burs
4002	Bullet - Jumbo 6 Burs

## POSTS

### ICELIGHT - LIGHT TRANSMITTING COMPOSITE POST

PART#	DESCRIPTION
90080-01	IceLight Intro Kit: 5 (1.0 mm) Posts, 5 (1.2 mm) Posts, 5 (1.4 mm) Posts, 5 (1.6 mm) Posts, 1 (1.0 mm) Reamer, 1 (1.2 mm) Reamer, 1 (1.4 mm) Reamer, 1 (1.6 mm) Reamer
90378-01	Refill (1.0 mm), 10 pack
90380-01	Refill (1.2 mm), 10 pack
90382-01	Refill (1.4 mm), 10 pack
90384-01	Refill (1.6 mm), 10 pack

### ICEPOST - COMPOSITE POST

PART#	DESCRIPTION
90079-01	IcePost Intro Kit : 5 (1.0 mm) Posts, 5 (1.2 mm) Posts, 5 (1.4 mm) Posts, 5 (1.6 mm) Posts, 1 (1.0 mm) Reamer, 1 (1.2 mm) Reamer, 1 (1.4 mm) Reamer, 1 (1.6 mm) Reamer
90082-01	Refill (1.0 mm), 10 pack

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## ICEPOST - COMPOSITE POST (CONTINUED)

PART#	DESCRIPTION
90083-01	Refill (1.2 mm), 10 pack
90088-01	Refill (1.4 mm), 10 pack
90089-01	Refill (1.6 mm), 10 pack

## REAMERS

PART#	DESCRIPTION
90325	Reamer (1.0 mm)
90327	Reamer (1.2 mm)
90329	Reamer (1.4 mm)
90331	Reamer (1.6 mm)

## PROVISIONAL COMPOSITES



## TURBOTEMP EZ

PART#	DESCRIPTION
90-00085	A-1 Refill (82 gm), 10 Mixing Tips
90-00086	A-2 Refill (82 gm), 10 Mixing Tips
90-00087	A-3/3.5 Refill (82 gm), 10 Mixing Tips
90-00088	B-1 Refill (82 gm), 10 Mixing Tips
90-00089	C-3 Refill (82 gm), 10 Mixing Tips
90-00090	Bleach White Refill (82 gm), 10 Mixing Tips

TURBOTEMP 2	
PART#	DESCRIPTION
90156	A-1 Refill (76 gm), 10 Mixing Tips
90157	A-2 Refill (76 gm), 10 Mixing Tips
92964	A-3 Refill (76 gm), 10 Mixing Tips
90158	A-3.5 Refill (76 gm), 10 Mixing Tips
90631	B-1 Refill (76 gm), 10 Mixing Tips
TURBOTEMP 3	
PART#	DESCRIPTION
93323	A-1 Refill (76 gm), 10 Mixing Tips
93324	A-2 Refill (76 gm), 10 Mixing Tips
93369	A-3 Refill (76 gm), 10 Mixing Tips
93325	A-3.5 Refill (76 gm), 10 Mixing Tips
93326	B-1 Refill (76 gm), 10 Mixing Tips
93591	Bleach Refill (76 gm), 10 Mixing Tips
TIPS AND ACCESSORIES	
PART#	DESCRIPTION
90162	Turbo Temp Mixing Tips, 10 pack
90163	Turbo Temp Mixing Tips, 50 pack
94746-15	TurboTemp EZ Mixing Tips, 15 pack
90176	Automix 4:1 Dispenser Gun
93844	Automix 10:1 Dispenser Gun
93199	Slider 10:1 Turbo Temp 3
TEMPORARY CEMENT	
PART#	DESCRIPTION
91087	Turbo Temp 2 Cement (10 ml)

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-TX
LOCATOR CHAIRSIDE
<b>MATERIALS &amp; EQUIPMENT</b>



TAKE-HOME KITS

PART#	DESCRIPTION
IVERISPEED	Iveri Speed: 1 LED light w/attachable tray, 4 Carbamide Peroxide (35%) Syringes (3.5 ml), Shade Chart





### TAKE-HOME REFILL SYRINGES

PART#	DESCRIPTION
4PKIVERI22	Take-Home Refill 22% Carbamide Peroxide - 4 Syringes (3.5 ml)
25PKIVERI22	Take-Home Refill 22% Carbamide Peroxide - 25 Syringes (3.5 ml)
4PKIVERI35	Take-Home Refill 35% Carbamide Peroxide - 4 Syringes (3.5 ml)
25PKIVERI35	Take-Home Refill 35% Carbamide Peroxide - 25 Syringes (3.5 ml)
4PKIVERI14	Take-Home Refill 14% Hydrogen Peroxide - 4 Syringes (3.5 ml)

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## AIR ABRASION



### PREPSTART UNITS AND ACCESSORIES

PART#	DESCRIPTION
91747	PrepStart H2O Hydro Abrasion Cavity Prep Unit
91913	PrepStart H2O Accessory Unit
200259-00	PrepStart Air Abrasion Unit
93383	Power Plus Air Booster
200628	PrepStart Handpiece Kit 1
91720	PrepStart H2O Handpiece & Nozzle (.019")
92944	PrepStart H2O Handpiece only

### PREPSTART NOZZLES

PART#	DESCRIPTION
186937-01	PrepStart Nozzle (.015"), 80°, 3 pack
186934-01	PrepStart Nozzle (.019"), 80°, 3 pack
91736	PrepStart H2O Nozzle (.019"), 80°
28515	PrepStart Tip (.015"), 45°
28528	PrepStart Tip (.015"), 80°
28520	PrepStart Tip (.015"), 90°
28521	PrepStart Tip (.015"), 120°
28519	PrepStart Tip (.019"), 45°
28527	PrepStart Tip (.019"), 80°

## PREPSTART NOZZLES (CONTINUED)

PART#	DESCRIPTION
28530	PrepStart Tip (.019"), 90°
28531	PrepStart Tip (.019"), 120°
28526	PrepStart Tip (.026"), 45°
28529	PrepStart Tip (.026"), 80°
28540	PrepStart Tip (.026"), 90°
28541	PrepStart Tip (.026"), 120°

## RETRACTORS

PART#	DESCRIPTION
89445	Standard Gingival Retractor
89448	Large Gingival Retractor

## SANDBLASTING



## MICROETCHER UNITS AND ACCESSORIES

PART#	DESCRIPTION
22005-01	MicroEtcher IIA Sandblaster
22000-03	MicroEtcher II Sandblaster
201936-00	MicroEtcher CD Kavo
201935-00	MicroEtcher CD Sirona
93825	MicroEtcher CD Bien Air
21000-03	MicroEtcher ERC
201502-00	Jar w/Grommet Lid (for MicroEtcher II)
25026	Jar w/Grommet Lid (for MicroEtcher ERC)

LOCATOR Root

LOCATOR  
Implants

LOCATOR  
Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL  
HYGIENE

LOCATOR R-TX

LOCATOR  
CHAIRSIDE

MATERIALS &  
EQUIPMENT

## MICROETCHER UNITS AND ACCESSORIES (CONTINUED)


PART#	DESCRIPTION
201503-01	MicroEtcher Filter, Jar Cap, 10 pack
22006	SandTrap 15 pack Disposable Dust Trap for MicroEtcher

## MICROETCHER NOZZLES



PART#	DESCRIPTION
201865-00	MicroEtcher IIA Nozzle Slim Profile (.048") , 60°
201802-00	MicroEtcher II Nozzle (.048") , 60°
186113	MicroEtcher II/IIA Nozzle Pit & Fissure (.032"), 60°
186002-01	MicroEtcher ERC Nozzle
186001	MicroEtcher II Standard Nozzle

## INSTALLATION KITS



PART#	DESCRIPTION
44000	Quick Disconnect Kit (for Operator)
44010	Quick Disconnect Second Station (for Operator)
44007	Quick Disconnect (for Lab)
44017	Quick Disconnect Second Station (for Lab)

## INSTALLATION KITS (CONTINUED)

PART#	DESCRIPTION
44020	Quick Disconnect Male
44075	MicroEtcher Handpiece Adapter 4-5 Hole
201431-00	MicroEtcher Handpiece Adapter Kavo Style

### SILJET

PART#	DESCRIPTION
93577	SilJet Kit: Siljet Powder (40 g), S-Bond (10 ml), E-Bond (10 ml), 40 applicator brushes, Tooth (2 g) & White (2 g) Accolade OP Mask, & 20 Tips (25 gauge)
93596	Powder Refill, 3 jars (40 g each)

### ABRASIVES



PART#	DESCRIPTION
80042A	Aluminum Oxide 27 Micron (1 lb.)
15301	Aluminum Oxide 50 Micron (1 lb.)
15201	Aluminum Oxide 90 Micron (1 lb.)
16101	Glass Beads 90 Micron (1 lb.)
14001	Orthoprophy SA-85 85 Micron (1 lb.)
17001	B+ Sodium Bicarb (MicroProphy) (1 lb.)
17011	Sodium Bicarb (other prophys) (1 lb.)

LOCATOR Root

LOCATOR Implants

LOCATOR Abutments

LOCATOR Bar

LOCATOR FIXED

REGAL HYGIENE

LOCATOR R-TX

LOCATOR CHAIRSIDE

MATERIALS & EQUIPMENT

## DUST CONFINEMENT



### DUST CABINETS

PART#	DESCRIPTION
91556	MacroCab Plus (120 v)
91010	MicroCab Plus (120 v)

### MICRO CAB REPLACEMENT PARTS

PART#	DESCRIPTION
49105	Replacement Window (MicroCab) Original
201359	10 Disposable Window Shields (MicroCab)
49055	Replacement Filter (MicroCab Plus)
59024	Handguard Kit (MicroCab)
91286	10 Disposable Window Shields (MicroCab) Original
93817	10 Disposable Window Shields (MicroCab) New

**MACROCAB REPLACEMENT PARTS**

PART#	DESCRIPTION
185508	Replacement Window (MacroCab)
201360	10 Disposable Window Shields (MacroCab)
116545	Replacement Filter (MacroCab)
59026	Handguard Kit (MacroCab)
93960	10 Disposable Window Shields (MacroCab Plus)
91557	Replacement Filter (MacroCab Plus)
91518	Handguard Kit (MicroCab Plus)

**DC-2 REPLACEMENT PARTS**

PART#	DESCRIPTION
49102	Replacement Window (DC-2)
49052	Replacement Filter (DC-2)

LOCATOR Root
LOCATOR Implants
LOCATOR Abutments
LOCATOR Bar
LOCATOR FIXED
REGAL HYGIENE
LOCATOR R-TX
LOCATOR CHAIRSIDE
<b>MATERIALS &amp; EQUIPMENT</b>

# RETURN POLICY AND WARRANTY

## ORDERING

Orders are accepted by internet, email, mail, fax, or phone. Regular business hours are 5:00am to 5:00pm PST (Pacific Standard Time). Domestic orders received by 1:00pm PST can be shipped by ground the same day. Orders received by 3:00pm PST can be shipped by express. Institutional orders require a purchase order number.

To place a new order, please choose one of the following ways:

1. Online: [www.zestdent.com](http://www.zestdent.com)
2. E-mail: [sales@zestdent.com](mailto:sales@zestdent.com)
3. Phone: 1-800-262-2310
4. Fax: 1-800-487-1357

## PRICING

Zest Dental Solutions makes every effort to maintain a competitive pricing structure. Pricing is subject to routine review and change without prior notice.

## TERMS AND BILLING

- E. Credit Card:** Payment at the time of order with VISA, AMERICAN EXPRESS, DISCOVER or MASTERCARD is necessary for customers not approved for open account billing or special financing offers.
- F. C.O.D.:** Delivery can also be made by C.O.D. freight collect.
- G. Open Account:** An open account can be established by completing a credit application and receiving approval by Zest Dental Solutions. The terms of orders purchased under open accounts are net 30 days.
- H. Export Orders:** Individual orders shipped outside the U.S. require payment in advance or a letter of credit.
- I. Past Due Accounts:** Past due balances will be subject to a 1.5% finance charge per month, amount equal to 18% per annum.

## SHIPPING CHARGES

Merchandise is shipped prepaid by Zest Dental Solutions with the cost added to the invoice, or freight collect with C.O.D. fee included in the case of C.O.D. shipments. Delivery options include Fed Ex ground, 3-day, 2-day, or overnight rush service.

## LIMITED WARRANTY

Zest Dental Solutions provides a limited warranty for its products, to the original purchaser, to be free from defects in workmanship and materials under normal use and service, for a period of one year from the date of purchase with the exception of LOCATOR and LOCATOR R-Tx abutments which carry a four year limited warranty (when purchased directly from Zest Dental Solutions) and a two year limited warranty when purchased from one of our partners. Zest Dental Solutions will, at its option, substitute the returned product that proves defective within the warranty period, with a similar product, free of charge.

Zest Dental Solutions continually strives to improve its products, and therefore, reserves the right to improve, modify, or discontinue products at any time without notice or incurring obligations. Purchaser assumes all risks and liability resulting from the use of Zest Dental Solutions Products, whether used separately or in combination with other products not of Zest Dental Solutions manufacture.

## RETURN POLICY

Please Observe the Following Guidelines:

- A. Authorization for returns must be received from Zest Dental Solutions prior to the return of any product. A Return Material Authorization (RMA) number will be provided for all returns.
- B. Shipping charges must be prepaid by the customer to accept a return shipment.
- C. Returned products are subject to a \$20 restocking fee for orders under \$400. For orders \$400 to \$1,000, a \$40 restocking fee will be applied. For orders exceeding \$1,000, please contact your Customer Service Representative for the restocking fee.
- D. Returned merchandise will be accepted within 90 days of purchase if product is in saleable condition (in its original unopened package and not marred by any added writing or over-labeling).
- E. Returns will not be accepted after 90 days of purchase.
- F. Non-returnable items include:
  7. Merchandise retained beyond expiration date noted on the package.
  8. Packages with broken seals or missing parts
  9. Used, damaged, or obsolete products will not be accepted for return.
- J. LOCATOR Implant returns will be accepted in a 1:1 exchange for a period of 6 months with no restocking fee.



# THE ZEST DENTAL SOLUTIONS® 10 YEAR WARRANTY



## THE ZEST DENTAL SOLUTIONS NARROW AND STANDARD RIDGE IMPLANT WARRANTY: LOCATOR® OVERDENTURE IMPLANT SYSTEM ("LODI").

Zest Dental Solutions ("Zest") is committed to providing quality products and dedicated to gathering feedback about its products. Zest actively collects and reviews the feedback of users of our products in compliance with regulatory reporting requirements and to better help us understand market expectations and validate our products' performance. The collection method Zest utilizes for such feedback is the Zest Product Experience Report ("PER") form. The PER form is to be completed from information provided by the attending clinician to share their Zest product experience.

Pursuant to the Zest implant warranty, Zest will replace covered LODI implants such qualifying covered implant is returned (as such terms are defined below). Upon a request for replacement under the Zest warranty, Zest will send a replacement product once Zest receives the returned product and completed PER form and confirms that the returned product is covered under the Zest warranty.

ONLY DIRECT ZEST CUSTOMERS MAY MAKE A WARRANTY CLAIM. In order to make a claim under the Zest warranty please return the items mentioned below in protective packaging and send these via a shipping method which enables the package to be tracked:

Printed copy of PER form completed from information provided by the attending clinician.

Explanted product(s) in sterile condition (NON-STERILE products will not qualify for replacement)

If a product failure or loss of integration has occurred, please additionally send relevant radiographs (these will not be returned unless specifically requested, please send copies)

Send shipment to: Zest Dental Solutions  
ATTN: Customer Service (US customers) or Wholesale Distribution  
(Distributors/OEM Partners)  
2875 Loker Ave East Carlsbad, CA 92010

If using the electronic form, please include the Tracking Number for the returned product package.

Download Zest NDI Product Experience Report (PDF) AT  
[WWW.ZESTDENT.COM](http://WWW.ZESTDENT.COM)

Zest Dental Solutions ("Zest") IMPLANT WARRANTY (Valid as of May 1, 2015)

**1. Warranty beneficiary and scope:** Zest Dental Solutions ("Zest") hereby warrants to the direct Zest customer purchasing the from Zest ("Customer") that the Zest LOCATOR Overdenture Implant ("LODI"), when implanted according to the respective Zest LODI Technique Manual and other written instructions provided by Zest by a clinician (the "User"), will be free from any loss or lack of integration, fracture or other structural failure for the period of 10 years ("Warranty Period") from the time of treatment by the User (collectively, the "Zest Warranty"). This warranty only applies to the Customer. Third parties, particularly patients, are not covered by the Zest warranty and have no rights hereunder. Customers' sole remedy and Zest's sole liability under this Zest warranty is the replacement of the LODI implant by Zest as set forth herein. The Zest implant warranty only covers the replacement of the LODI implant and not any associated costs or expenses, including, but not limited to, chair time, laboratory fees and any other associated treatment.

**2. Zest Warranty conditions:** In the event that any request for warranty service for an implant is made by a Customer under the Zest warranty during the Warranty Period, Zest will replace such implant with the same or substantially equivalent product subject to the terms and conditions herein. The replacement product will be sent upon receipt of the returned product and the completed PER form and confirmation that the product is covered under the Zest warranty. To qualify for coverage under the Zest warranty, the claim must be made within the applicable Warranty Period by a Customer and all conditions below must be met. Once these requirements are satisfied, Zest will send the replacement part(s). The following REQUIRED conditions must be met and documented in order for coverage of a returned implant under

this Zest implant warranty:

1. The LODI was used exclusively with all components, connections, attachments and other technology provided by Zest and not in combination with any other manufacturer's products or technology;
2. The LODI is returned in sterile condition (or disinfected if delivered as such);
3. The LODI was implanted by a User and inspected and maintained in full compliance with the respective Zest LODI Technique Manual valid at the time of treatment and all other Zest written instructions as well as recognized dental procedures, during and after the treatment;
4. The patient had good oral hygiene which was monitored and documented by the User
5. The LODI was not subjected to damage caused by misuse, misapplication, accident, trauma or any other damage caused by external factors or the User, the patient or a third party;
6. A completed and signed PER was completed and submitted no later than 10 days after the complaint is made. Customer is responsible to ensure that the PER is completed and submitted from information provided by the applicable User. Details of the incident are imperative to determine whether vigilance reporting is required to regulatory authorities.

**3. Limits and limitations:** This Zest implant warranty is the only guarantee provided by Zest and shall apply in addition to any warranty rights conferred under any written sales agreement executed by Zest. The User remains free to claim rights against his supplier. EXCEPT AS SET FORTH HEREIN AND IN A WRITTEN SALES AGREEMENT EXECUTED BY ZEST, ZEST HEREBY DISCLAIMS ANY OTHER WARRANTIES, EXPRESS OR IMPLIED, WITH RESPECT TO LODI/SNDI OR ANY OTHER ZEST PRODUCTS, SERVICES OR INFORMATION, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTY OF FITNESS FOR PURPOSE OR MERCHANTABILITY, OR NON-INFRINGEMENT OF THIRD PARTY RIGHTS, OR ANY WARRANTY BY COURSE OF PERFORMANCE OR OTHERWISE.

ZEST SHALL NOT BE LIABLE TO THE CUSTOMER, THE USER, THE PATIENT, OR ANY THIRD PARTY, FOR LOST EARNINGS OR PROFITS, DIRECT OR INDIRECT DAMAGES AS WELL AS COLLATERAL, INCIDENTAL, PUNITIVE OR CONSEQUENTIAL DAMAGES, DIRECTLY OR INDIRECTLY RELATED TO LODI OR ANY OTHER ZEST PRODUCTS, SERVICES OR INFORMATION.

4. Zest Warranty territory: This Zest implant warranty applies worldwide to LODI implants sold by Zest, a Zest affiliated company or an official distributor of Zest on or after the validity date stated above.
5. Modification or termination: Zest may modify or terminate this Zest implant warranty at any time in whole or in part. Changes to, or the termination of the Zest implant warranty, will not affect the warranty given for LODI installed prior to the date of the change or termination.

### LOCATOR OVERDENTURE IMPLANT (LODI) PRODUCT EXCHANGE POLICY

Zest Dental Solutions understands that customers may need to adjust inventories of LOCATOR Overdenture Implant System (LODI) Products in order to achieve the correct mix of sizes to treat their patients. Zest Dental Solutions will waive normal restocking fees (1:1 or greater) for exchanges during the first six months following the original purchase. Packaging cannot be written on or in any way adulterated. Shipping will still be the responsibility of the customer requesting the exchange.





# CONTACT US

1.800.262.2310

SALES@ZESTDENT.COM

LEARN MORE AT  
[WWW.ZESTDENT.COM](http://WWW.ZESTDENT.COM)



GLOBAL HEADQUARTERS  
2875 LOKER AVE EAST  
CARLSBAD, CA 92010  
USA